


L-column ODS (3 μ m)

L-column ODS(3 μ m) is also produced using a preeminent deactivation processing technology just like L-column ODS(5 μ m). The preeminent deactivation is performed by high-temperature gaseous-phase silylation (Super endcapping patent No.5,134,110 U.S.A., No.2611545 JAPAN). It deactivates the residual silanol on the silica gel surface nearly completely.

Theoretical plate number with L-column ODS(3 μ m) can be obtained at 1.5 times (150,000 plates/meter) the rate with which they are obtained with the L-column ODS(5 μ m), which is the highest that has been obtained so far for commercial columns. L-column ODS(3 μ m) is suitable for high-throughput analysis because its optimal flow rate is twice higher than that of 5 μ m column.



Specifications of Silica particles

Average particle size: 3 μ m
 Average pore size: 12nm
 Surface area: 340m²/g
 Pore volume: 1.1mL/g
 Carbon content: 17%


■ Anti-depressants' analysis

Anti-depressants are basic compounds. In HPLC analysis, they cause tailing of the peak. We attempted simultaneous analysis of typical anti-depressants with the L-column ODS(3 μ m).

There was no tailing of the peak, and adsorption occurred without an ion-pair reagent.

Since L-column ODS(3 μ m) yields theoretical plate number at 1.5 times the rate at which they are obtained with L-column ODS(5 μ m), the analysis time is shortened and a column length of two-thirds yields same separation efficiency.

Therefore it can be said that this is the most suitable column for highly sensitive separation of basic compounds and high-speed analysis.





[Analytical conditions]

Column : a) L-column ODS (5 μ m) 4.6 x 150mm
 b) L-column ODS (3 μ m) 4.6 x 100mm
 Mobile Phase : Methanol/50mM CH₃COONH₄=66/34
 Flow-Rate : 1mL/min
 Temp. : 30°C
 Detection : UV 254 nm

■ Relationship between theoretical plates and the flow-rate

The optimal flow-rate of a 2.1 mm i.d. semi-micro column packed 5 μ m ODS is 0.2 mL/min. Since L-column ODS (3 μ m) can yield a maximum theoretical plate number rate of 0.25 – 0.4 mL/min, it is probably the best column most suitable for shortening of time and high-throughput analysis. Furthermore, it is also effective for LC/MS analysis using a short column (50 mm and 100mm), because its theoretical plate number is better than those of columns produced by other companies.



■ High-speed analysis of "cold" medicines

We analyzed six ingredients in a commercially available "cold" medicine. When a 4.6x150mm L-column ODS (5 μ m) was changed to a 4.6x50mm L-column ODS (3 μ m), the six ingredients were separated completely in less than 2 minutes when the flow-rate was set at 2 mL/min.

[Analytical conditions]

Column : a) L-column ODS (5 μ m) 4.6 x 150mm
 b) L-column ODS (3 μ m) 4.6 x 50mm
 Mobile Phase : Methanol/20mM H₃PO₄ = 40/60
 Flow-Rate : a) 1mL/min
 : b) 2mL/min
 Temp. : 40°C
 Detection : UV 254 nm

1. Acetaminophen
2. Caffeine
3. Narcotine
4. Chlorpheniramine Maleate
5. Ethenzamide
6. Bucetin

