

HPLC column

# Sunniest



**Sunniest C18-HT 2  $\mu$ m**

# Sunniest C18-HT, 2 $\mu$ m


## Features

- Low back pressure and high efficiency by precisely classified particle
- High pressure packing (10,000 psi) using hard silica gels with high pressure resistant • leads long column life without any void.
- Unique bonding technique for Sunniest (patent pending)
- The most suitable inner surface of column by special grinding

**Sunniest technique  
to 2 $\mu$ m particle!**


## • Narrow Particle Distribution and Low Back Pressure

Measured by Coulter Counter method




Conventional 2  $\mu$ m silica gel particle was classified again. 20% volume was cut off from both sides respectively. Consequently column back pressure reduced more than 15%. Our 2  $\mu$ m silica gel particle shows a half pressure to compare with the other sub-2  $\mu$ m silica gel particle.

Column pressure using methanol or acetonitrile and water



Column: Sunniest C18-HT, 2  $\mu$ m 50 x 2.1 mm  
Mobile phase: CH<sub>3</sub>OH/H<sub>2</sub>O, CH<sub>3</sub>CN/H<sub>2</sub>O  
Flow rate: 0.5 mL/min  
Temperature: 40  $^{\circ}$ C

## • An Unique Modification (Patent Pending as Sunniest series)



Column: Sunniest C18-HT, 2  $\mu$ m 50 x 2.0 mm

Mobile phase:

**A)** CH<sub>3</sub>OH/20mM Phosphate buffer pH7.5 = 80/20


**B)** CH<sub>3</sub>CN/20mM Phosphate buffer pH7.0 = 60/40

Flow rate: 0.4 mL/min

Pressure: **A)** 19.5 MPa, **B)** 13.5 MPa


Temperature: 40  $^{\circ}$ C

Sample: 1 = Uracil, 2 = Propranolol,



It is difficult to end-cap on sub 2  $\mu$ m or 2  $\mu$ m silica gel particle as well as 3  $\mu$ m or 5  $\mu$ m silica gel particle. Most sub 2  $\mu$ m or 2  $\mu$ m C18 columns show smaller plate number and higher tailing factor than Sunniest C18-HT. Sunniest C18-HT 2  $\mu$ m shows good peak shape for amitriptyline under not only methanol/phosphate buffer mobile phase but also acetonitrile/phosphate buffer mobile.

## • Separation of Analgesics



Mobile phase: CH<sub>3</sub>CN/0.1% Formic acid = 20/80  
Flow rate: 1.0 mL/min for 150 x 4.6 mm


**0.6 mL/min for 50 x 2.1 mm**

Temperature: 40 °C

Detection: UV@230 nm


Sample:

1=Acetaminophen, 2=Antipyrine, 3=Aspirin, 4=Ethenzamide



2  $\mu$ m particle allows to reduce retention time because high efficiency is kept under high flow rate conditions. As shown the above chromatograms, analytical time reduced 1/8 without sacrifices of separation by using 2  $\mu$ m, 50 x 2.1 mm column instead of 5  $\mu$ m 150 x 4.6 mm column.

## • Comparison of Plate Number



Mobile phase: CH<sub>3</sub>CN/H<sub>2</sub>O = 60/40

Flow rate: 0.6 mL/min for 2.1 x 30 mm and 2.1 x 50 mm

1.0 mL/min for 3.0 x 50 mm

0.4 mL/min for 2.1 x 75 mm and 2.1 x 100 mm

Temperature: 40 °C

Detection: UV@250 nm

Sample: 1=Uracil,

2=Toluene,

3=Acenaphthene,

4=Butylbenzene

## • Characteristics of Sunniest C18-HT, 2 $\mu$ m

Packings	Silica gel support			C18			
	Particle size ( $\mu$ m)	Pore diameter (nm)	Specific surface area (m <sup>2</sup> /g)	Carbon content (%)	Bonded phase	Maximum operating pressure	Available pH range
Sunniest C18-HT	2.0 (Coulter counter)	10	340	16	C18	70 MPa or 10,000 psi	1.5 - 10

It is very important for 2 mm particle to have a capacity to resist pressure because of high column back pressure. The larger a pore volume of silica gel, the weaker a capacity to resist pressure. The silica gel with 0.85 ml/g of pore volume is used for Sunniest C18-HT, 2 mm, so that it have a high capacity to resist pressure and a high operating pressure.




## • Surface Roughness on Inner Surface of Column

### Parameter of surface roughness

Ra: Average roughness from center line


Rz: Roughness calculated from 10 points average (5 points of maximum and 5 points of minimum)



Schematic diagram of surface roughness

	G company	Y company	S1 company	S2 company	W company	ChromaNik Technologies
Ra	0.34 $\mu\text{m}$	0.32 $\mu\text{m}$	0.37 $\mu\text{m}$	0.03 $\mu\text{m}$	0.20 $\mu\text{m}$	0.01 $\mu\text{m}$
Rz	1.88 $\mu\text{m}$	1.62 $\mu\text{m}$	1.91 $\mu\text{m}$	0.19 $\mu\text{m}$	0.90 $\mu\text{m}$	0.10 $\mu\text{m}$

It is considered that surface roughness affects column performance. Surface asperity of ChromaNik Technologies column is 1/30 to 1/20 to compare with that of GL Sciences, YMC, Shimadzu and Waters columns. ChromaNik Technologies provides a column with very smooth surface which is the most suitable for 2  $\mu\text{m}$  particle packing.



## • Ordering Information of Sunniest C18-HT, 2 $\mu\text{m}$

Length (mm)	30	50	75	100
Inner diameter (mm)	Cat. No.	Cat. No.	Cat. No.	Cat. No.
2.1	EB1931	EB1941	EB1951	EB1961
3.0	EB1331	EB1341	EB1351	EB1361

### Manufacturer

### ChromaNik Technologies Inc.

6-3-1 Namiyoke, Minato-ku, Osaka, 552-0001 Japan

TEL: +81-6-6581-0885 FAX: +81-6-6581-0890

E-mail: [info@chromanik.co.jp](mailto:info@chromanik.co.jp)

URL: <http://chromanik.co.jp>