


Heraeus



Vacuum UV Light Sources
for Laboratory and Production Applications

Vacuum UV Light Sources

30 Watt Types

Heraeus Noblelight's vacuum UV light sources encompass the widest range of 30 and 200 Watt deuterium lamps on the market. While the 30 Watt lamps are used as calibration sources in the far UV (V-series) and for analytical purposes, the 200 Watt water-cooled high power deuterium lamp is a source of high photon energy deep ultraviolet radiation. It delivers a radiant intensity 4 to 5 times higher than standard 30 Watt deuterium lamps.

Vacuum UV Light sources, 30 Watt type - V series

The V series lamps consist of five high stability deuterium arc sources for use as a calibration source in the far ultraviolet. The 10 Volt heater, 1 mm arc aperture lamps are designed and built with a rugged construction with magnesium fluoride windows. Output extends to the Lyman-alpha line and below.

Application: VUV lamp for all applications that require high stability output down to 115 nm.

General Construction

There are two overall lengths options: The longer version is built for the highest light output stability. The light is prevented from being reflected from the sides of the mounting tubes by multiple baffles.

Both versions of the lamp are available either with a synthetic silica window for emission from 160 nm upwards or with a MgF₂ (magnesium fluoride) window that allows emission from 115 nm upwards.



V03 Light Source

Specifications 30 Watt VUV Lamp Models

| Lamp Type | V01 | V02 | V03 | V04 | V05* |
|----------------------------|------------------|-------------|------------------|-------------|------------------|
| Part number | 80017697 | 80017698 | 80017699 | 80017700 | 80020174 |
| Spectral range | 115-400nm | 160-400nm | 115-400nm | 160-400nm | 115-400nm |
| Window material | MgF ₂ | Spectrosil | MgF ₂ | Spectrosil | MgF ₂ |
| Aperture diameter | 1.0 mm | 1.0 mm | 1.0 mm | 1.0 mm | 1.0 mm |
| Arc height from window | 102 mm | 102 mm | 72 mm | 72 mm | 102 mm |
| Overall length incl. base | 146-150 mm | 146-150 mm | 116-120 mm | 116-120 mm | 140-144 mm |
| Dimensional outline | see page xy | | | | |
| Electrical data | | | | | |
| Heater voltage start | 10Vdc | 10Vdc | 10Vdc | 10Vdc | 10Vdc |
| Heater voltage operational | 0/6 V | 0/6 V | 0/6 V | 0/6 V | 0/6 V |
| Heater current | 1.0Adc max. | 1.0Adc max. | 1.0Adc max. | 1.0Adc max. | 1.0Adc max. |
| Anode voltage | 60-80Vdc | 60-80Vdc | 60-80Vdc | 60-80Vdc | 60-80Vdc |
| Anode current | 300mAdc | 300mAdc | 300mAdc | 300mAdc | 300mAdc |
| Light output | | | | | |
| Noise | ≤ 0.2% p-p | ≤ 0.2% p-p | ≤ 0.2% p-p | ≤ 0.2% p-p | ≤ 0.2% p-p |
| Drift | ± 0.5%/h | ± 0.5%/h | ± 0.5%/h | ± 0.5%/h | ± 0.5%/h |
| Lifetime (warranty) | ≥ 500h | ≥ 500h | ≥ 500h | ≥ 500h | ≥ 500h |

*Type V05 has less base and is specifically designed for enclosed vacuum operation

High Power VUV Light Sources

200 Watt Types



The D 200 high power VUV light sources are water-cooled deuterium discharge lamps of 200 Watt electrical power. They deliver a radiant intensity 4-5 times higher than standard 30 Watt deuterium lamps. To remove the excess heat of the lamp the bulb is enclosed in a water cooled metal fixture engineered to take most of the heat produced and transfer it through water to a designated cooling reservoir. The D 200 F lamp has a synthetic silica window for an efficient emission of UV radiation between 160 nm and 400 nm, while the D 200 VUV uses an MgF₂ window to allow vacuum UV radiation down

to 115 nm. This lamp is provided with a vacuum flange DN 50 KF (or optional DN 63 CF-flange) for mounting on a vacuum chamber.

Applications

- Optical research with UV radiation
- UV spectroscopy at high radiation density
- Removal of electrostatic charge from semiconductor wafer
- Photochemical processes
- UV absorption measurements
- Fluorescence excitation


Specifications D 200 Lamp Types

| Lamp Type | D 200 F | D 200 F-HV | D 200 VUV |
|--|--------------------------------|------------------|--------------------------|
| Part number, with cooling jacket | 56001671 | 45006278 | 45006010 |
| Spectral range | 160-400 nm | 160-400 nm | 115-400 nm |
| Window material | Synthetic silica | Synthetic silica | MgF ₂ |
| Vacuum flange | no | DN 50 KF | DN 50 KF (optional CF63) |
| Aperture diameter | 1.0 mm | 1.0 mm | 1.0 mm |
| Arc height from window | | | |
| Overall length | 158 mm | 243 mm | 243 mm |
| Outer diameter | 54 mm | 54 mm | 54 mm |
| Cooling/flow rate | Water cooling / ≥0.5 l/min. | | |
| Cooling water pressure | ≤ 0.3 bar | ≤ 0.3 bar | ≤ 0.3 bar |
| Connector | Amphenol-Tuchel 3109/1 | | |
| Weight | 0.65 kg | 0.9 kg | 0.9 kg |
| Dimensional outline | See page 6 | | |
| Replacement lamp part no. without cooling jacket | 56001669 | 45006286 | 45006127 |
| Electrical data | | | |
| Heater voltage-start | 6.0Vdc | 6.0Vdc | 6.0Vdc |
| Heater current-start | 4.5Adc | 4.5Adc | 4.5Adc |
| Warm-up time | 30s | 30s | 30s |
| Heater voltage operational | 3.0Vdc | 3.0Vdc | 3.0Vdc |
| Heater current operational | 2.5Adc | 2.5Adc | 2.5Adc |
| Strike voltage | ≥ 500 V | ≥ 500 V | ≥ 500 V |
| Anode voltage | 110 +10/-15Vdc | 110 +10/-15Vdc | 110 +10/-15Vdc |
| Anode current | adjustable 0.9 to 1.8Adc | | |
| Light output | | | |
| Noise @ 250nm | ≤ 0.1 % p-p | ≤ 0.1 % p-p | ≤ 0.1 % p-p |
| Drift @ 250nm | ± 0.5 %/h | ± 0.5 %/h | ± 0.5 %/h |
| Lifetime @ 280nm | ≥ 1000 h at 0.9 A lamp current | | |


Vacuum Flange DN 63 CF for D 200 VUV

To make use of the UV radiation below 180nm the lamp D 200 VUV is often used on high vacuum chambers. For this purpose Heraeus provides the D 200 VUV lamp with a DN 50 KF flange as a standard flange. In case the vacuum chamber is equipped with CF flanges Heraeus offers a DN 63 CF adapter flange. Other CF-flange dimensions can be offered on request. (dimensional outline see last page).

Spectral Energy Distribution in VUV range



Spectral Energy Distribution in UV range



Power Supply for D200 VUV Lamps


The power supply PSD 200 is designed to operate all above mentioned D 200 Lamp types. It works automatically. If required the lamp output can be varied by adjusting the anode current on the front panel between 0.9 to 1.8Adc. The adjusted current can be read on a LCD display. Also shown on a display is the elapsed operating time of the lamp.




Specifications PSD 200 Power Supply

| Power Supply | PSD 200 |
|---|--|
| Anode voltage | 90-145Vdc |
| Anode current | adjustable between 0.9-1.8Adc |
| Strike voltage | 600 V |
| Heater voltage start | 6.0 Vdc / 4.5 Adc |
| Heater voltage operational | 3.0 Vdc / 2.5 Adc |
| Mains input (AC) | 100/115/230 Vac, 50/60 Hz |
| Power consumption | 250 VA |
| Operating ambient temperature | 0°C to 40°C |
| Cooling method | closed loop water circulator, flow ≥ 0.5 l/min |
| Protection of lamp in case of malfunction | by flow switch if less cooling water, water circulator with flow switch to be provided by the customer |
| Dimensions (W x H x D) | 21 x 17.5 x 30 cm |
| Weight | 12 kg |
| Part Numbers | |
| Mains input: 230/115 V, 50/60Hz | 45006028 |
| Mains input: 100V, 50/60 Hz | 45006277 |


Dimensional Outlines




V01+V02




V03+V04





V05



D200F-HV, D200VUV



D200F



Flange DN 63 CF for D200 VUV

Europe, Middle East, Africa, Rest of World*

Heraeus Noblelight GmbH

Heraeusstrasse 12-14

63450 Hanau, Germany

Phone +49 6181 35 5086

Fax +49 6181 35 7970

hng-analyticalamps@heraeus.com

www.heraeus-noblelight.com

America*

Heraeus Noblelight LLC

1520C Broadmoor Blvd.

Buford 30518, GA, USA

Phone +1 678 835 5681

Fax +1 678 835 5766

sales.hni@heraeus.com

www.heraeus-noblelight.com

Asia-Pacific, Oceania*

Heraeus Noblelight (Shenyang) Ltd.

Room 502, 5F, 16th building

No. 99 Tianzhou Road

200233 Shanghai, PR China

Phone +86 21 5445 2255

Fax +86 21 5445 2410

info.hns@heraeus.com

www.heraeus-noblelight.cn

*For local contacts please visit also our website <http://www.heraeus-noblelight.com/en/contact/worldmap.aspx>