



Analysis of Pesticides, Mycotoxins, and Cannabinoids in Cannabis Gummies

By Nathaly Reyes, PhD

Abstract

The state of California demands the analysis of pesticides, mycotoxins, and cannabinoids in all types of cannabis-derived goods [1]. For that reason, reliable workflows for the easy determination of these analytes in diverse matrices are highly desired. Gummies are a very popular cannabis edible, and their composition makes them a highly complex matrix. In this work, we describe a complete workflow for the analysis of California pesticides, mycotoxins, and cannabinoids in gummies using a single extraction procedure.

Introduction

Edibles infused with cannabis or cannabidiol (CBD) are growing in popularity with cannabis consumers. In terms of testing, potency analysis of edibles is mandatory in all states where lab testing is a requisite. In addition, in the state of California, determination of contaminants, such as pesticides and mycotoxins, in edibles and other cannabis-derived goods is also a requirement. Due to the broad variety of edibles available in the market, different analytical strategies should be pursued in order to obtain reliable analytical data for all analytes of interest in all matrices. Among the various types of cannabis edibles, gummies present unique challenges for cannabis testing labs.

Gummies are a type of sticky matrix typically made of ingredients such as sugar, starch, pectin, and gelatin. Due to their composition, sample preparation strategies for potency testing usually include a sample solubilization step in solvents like water or DMSO or, alternatively, large volumes of solvents like methanol may be used. Following sample homogenization, cannabinoids are quantified via HPLC-UV by injecting diluted extract either directly or following a salting-out step using QuEChERS salts. As for the analysis of contaminants like pesticides and mycotoxins in gummy matrix, publicly available information is quite scarce, with only one report showing semiquantitative data for 35 pesticides using QuEChERS [2].

In this work, we provide a robust workflow for the quantitative determination of the California list of pesticides, mycotoxins, and cannabinoids in gummy samples. The optimized sample preparation methodology involves sample solubilization followed by an extraction step using acidified acetonitrile, and a salting-out step using EN QuEChERS salts. For the analysis of cannabinoids and LC-amenable contaminants, a simple dilution was conducted prior to injection. Whereas, for the GC-amenable pesticides, the use of a dSPE sorbent mix including primary/secondary amine (PSA), graphitized carbon black (GCB), and magnesium sulfate was necessary before analysis. Overall, the proposed workflow for the analysis of pesticides, mycotoxins, and cannabinoids in cannabis gummies provides satisfactory results in terms of linearity, accuracy, precision, and limits of quantitation (LOQs).


Experimental

Sample Preparation

After preliminary testing of various experimental parameters, the following workflow was implemented (Figure 1). Chopped gummies (1 g) were weighed into a 50 mL tube (cat.# 25846), and 5 mL of water was added. Then, the samples were vortexed vigorously until the gummy pieces were solubilized. The sample solution was then fortified with target analytes and/or internal standards (internal standards were used only for contaminant analysis, not for potency testing) as appropriate and vortexed again for 30 seconds. Next, 5 mL of acetonitrile acidified with acetic acid (1% v:v) was added, and the samples were vortexed another 30 seconds. A pouch of European EN 15662 Q-sep QuEChERS extraction salts (cat.# 25849) was added, and the sample was vortexed again and then centrifuged for 5 minutes.

For potency testing, 100 μ L of the sample extract was mixed with 900 μ L of 25:75 water:acetonitrile, and 2 μ L were analyzed by LC-UV. For the LC-amenable contaminants, 750 μ L of the sample extract was mixed with 250 μ L of water, and 2 μ L were analyzed by LC-MS/MS. For the GC-amenable contaminants, 1.9 mL of extract was transferred to a Q-sep QuEChERS dSPE tube containing pre-weighed PSA, GCB, and magnesium sulfate sorbents (cat.# 26217). After vortexing and centrifuging, 500 μ L of extract was mixed with 500 μ L of 1% acetic acid in acetonitrile, and 1 μ L was analyzed by GC-MS/MS.

Figure 1: Sample preparation procedure for the analysis of pesticides, mycotoxins, and cannabinoids in cannabis gummies.



Quantitation

For the analysis of pesticides and mycotoxins, calibration solutions were prepared by spiking analytes and internal standards in aliquots of extract obtained from blank gummy samples (a pooled extract was obtained by mixing extracts from various blank gummy samples). Table I shows the volume of target analyte spiking solution and internal standard mix solution that was added to each aliquot (final calibration solution volume = 3 mL). To construct calibration curves for the GC-amenable pesticides, 1.9 mL of the 3 mL calibration solutions were subjected to dSPE cleanup as described in the sample preparation section. Method accuracy and precision were evaluated by spiking homogenized gummies (after adding water and vortexing) at 10, 50, 100, and 500 ng/g in quadruplicate (Table II), and performing the full sample preparation workflow described in the sample preparation section.

For the analysis of cannabinoids, calibration solutions were prepared at 2, 5, 10, 20, 50, 100, and 200 ppm in 75:25 acetonitrile:water. Recovery of the cannabinoids was assessed by spiking gummy samples (1 g) solubilized in water at 0.2 and 0.5 mg/g (n=2), and then extracting them as previously described.

Table I: Preparation of calibrators for pesticides and mycotoxins analysis using aliquots of extracts collected from blank gummy samples (final volume of each calibration solution = 3 mL).

| Desired Analyte Conc. in Matrix (ng/g) | Analyte Conc. in Final Extract Assuming 100% Recovery from Matrix (ng/mL) | μL of Target Analyte Solution Spiked into Blank Extract | Analyte Conc. in Spiking Solution (ng/mL) | μL of 5000 ng/mL Internal Standard Mix Added to Blank Extract |
|--|---|--|---|--|
| 5 | 1 | 30 | 100 | 24 |
| 20 | 4 | 12 | 1000 | 24 |
| 50 | 10 | 30 | 1000 | 24 |
| 75 | 15 | 45 | 1000 | 24 |
| 150 | 30 | 90 | 1000 | 24 |
| 200 | 40 | 24 | 5000 | 24 |
| 400 | 80 | 48 | 5000 | 24 |
| 700 | 140 | 84 | 5000 | 24 |

Table II: Fortification of mycotoxins and pesticides in gummies (1 g of sample homogenized with 5 mL of water) at different concentration levels.

| Conc. in Matrix (ng/g) | μL of Target Analyte Solution Spiked into Homogenized Samples in Water | Analyte Conc. in Spiking Solution (ng/mL) | μL of 5000 ng/mL Internal Standards Mix Added to Homogenized Samples in Water |
|------------------------|---|---|--|
| 10 | 10 | 1000 | 40 |
| 50 | 50 | 1000 | 40 |
| 100 | 100 | 1000 | 40 |
| 500 | 100 | 5000 | 40 |

Instrument Conditions

Instrumentation and conditions for the analysis of pesticides, mycotoxins, and cannabinoids in cannabis gummies are presented in Tables III, IV, and V. Ion transitions for LC-amenable and GC-amenable contaminants are presented in Tables VI and VII, respectively. Cannabinoid retention times are shown in Table VIII.

Table III: LC-MS/MS Conditions (Pesticides and Mycotoxins).

| | | | | |
|-------------------|---|-----|------------|-----|
| Column | Raptor ARC-18 2.7 µm, 100 mm x 2.1 mm (cat.# 9314A12) | | | |
| Guard Column | Raptor ARC-18 EXP guard column cartridge 2.7 µm, 5 x 2.1 mm (cat.# 9314A0252) | | | |
| Mobile Phase A | Water, 2 mM ammonium formate, 0.1% formic acid | | | |
| Mobile Phase B | Methanol, 2 mM ammonium formate, 0.1% formic acid | | | |
| Time Program | Time (min) | %B | Time (min) | %B |
| | 0 | 5 | 10.5 | 100 |
| | 1.5 | 65 | 10.6 | 5 |
| | 8.5 | 95 | 12.0 | 5 |
| | 9.5 | 100 | — | — |
| Flow | 0.5 mL/min | | | |
| Column Temp. | 40 °C | | | |
| Autosampler Temp. | 10 °C | | | |
| Inj. Volume | 2 µL | | | |
| Instrument | Shimadzu LCMS-8060 | | | |

Table IV: GC-MS/MS Conditions (Pesticides).

| | |
|---------------------|--|
| Instrument | Thermo Trace 1310-TSQ 8000 |
| Column | Rxi-5ms, 30 m x 0.25 mm ID x 0.25 µm, (cat.# 13423) |
| Injection Mode | Splitless |
| Inj. Vol. | 1 µL |
| Liner | Topaz 4.0 mm ID single taper inlet liner w/ wool (cat.# 23447) |
| Inj. Temp. | 250 °C |
| Split Flow | 14.0 mL/min |
| Splitless Time | 0.50 min |
| Purge Flow | 5 mL/min |
| Oven | 90 °C (hold 1 min) to 310 °C at 25 °C/min C (hold 10 min) |
| Carrier Gas | He, constant flow |
| Flow Rate | 1.40 mL/min |
| Detector | MS/MS |
| Method Type | Acquisition - timed |
| Ionization Mode | EI |
| Transfer Line Temp. | 290 °C |
| Source Temp. | 330 °C |

Table V: HPLC-UV Conditions (Potency Using a Solvent Savings Method [3]).

| | |
|-------------------|---|
| Instrument | Waters ACQUITY |
| Column | Raptor ARC-18 2.7 µm, 150 mm x 2.1 mm (cat.# 9314A62) |
| Guard column | Raptor ARC-18 EXP guard column cartridge 2.7 µm, 5 x 2.1 mm (cat.# 9314A0252) |
| Inj. Vol. | 2 µL |
| Mobile phase A | Water, 5 mM ammonium formate, 0.1% formic acid |
| Mobile phase B | Acetonitrile, 0.1% formic acid |
| Gradient | Isocratic, 75% B |
| Flow | 0.4 mL/min |
| Column Temp. | 30 °C |
| Autosampler Temp. | 10 °C |
| Wavelength | 228 nm |

Table VI: LC-MS/MS Transitions.

| Name | Retention Time (min) | Precursor Ion | Product Ion 1 | Product Ion 2 |
|---------------------|----------------------|---------------|---------------|---------------|
| Daminozide-D6 | 0.7 | 167.0 | 149.3 | 49.3 |
| Daminozide | 0.7 | 161.1 | 44.1 | 143.2 |
| Acephate | 1.7 | 184.0 | 143.1 | 95.1 |
| Oxamyl | 2.0 | 237.1 | 72.1 | 90.1 |
| Flonicamid | 2.1 | 230.1 | 203.1 | 174.1 |
| Methomyl | 2.1 | 163.1 | 88.1 | 106.1 |
| Thiamethoxam | 2.1 | 292.0 | 211.1 | 181.1 |
| Imidacloprid | 2.3 | 256.1 | 209.1 | 175.1 |
| Mevinphos | 2.4 | 225.1 | 127.1 | 193.2 |
| Acetamiprid | 2.4 | 223.0 | 126.1 | 56.1 |
| Dimethoate-D6 | 2.4 | 236.1 | 205.1 | - |
| Dimethoate | 2.4 | 230.0 | 199.1 | 125.1 |
| Thiacloprid | 2.5 | 253.0 | 126.0 | 90.1 |
| Aflatoxin G2 | 2.5 | 331.2 | 189.3 | 115.2 |
| Aflatoxin G1 | 2.5 | 329.2 | 243.2 | 215.3 |
| Aldicarb | 2.6 | 116.0 | 89.2 | 70.2 |
| Aflatoxin B2 | 2.6 | 315.3 | 287.2 | 243.3 |
| Dichlorvos | 2.7 | 220.9 | 109.1 | 79.2 |
| Dichlorvos-D6 | 2.7 | 227.0 | 115.1 | - |
| Aflatoxin B1 | 2.7 | 313.2 | 241.2 | 128.2 |
| Imazalil | 2.7 | 297.0 | 159.0 | 201.0 |
| Carbofuran | 2.7 | 222.1 | 123.1 | 165.2 |
| Propoxur | 2.7 | 210.1 | 111.1 | 93.1 |
| Carbaryl-D7 | 2.8 | 209.2 | 152.2 | - |
| Carbaryl | 2.8 | 202.1 | 145.1 | 127.1 |
| Diuron-D6 | 3.0 | 239.1 | 78.2 | - |
| Atrazine-D5 | 3.0 | 221.2 | 179.1 | - |
| Naled | 3.1 | 397.8 | 127.1 | 109.1 |
| Metalaxyl | 3.1 | 280.2 | 220.2 | 192.2 |
| Spiroxamine | 3.1 | 298.3 | 144.2 | 100.2 |
| Chlorantraniliprole | 3.2 | 483.9 | 452.9 | 285.9 |
| Phosmet | 3.2 | 318.0 | 160.1 | 77.2 |
| Azoxystrobin | 3.3 | 404.0 | 372.1 | 344.1 |
| Linuron-D6 | 3.3 | 255.1 | 160.1 | - |
| Fludioxonil* | 3.4 | 247.0 | 180.0 | 126.0 |
| Methiocarb | 3.4 | 226.1 | 169.1 | 121.1 |
| Dimethomorph | 3.5 | 388.2 | 301.2 | 165.3 |
| Boscalid | 3.5 | 342.9 | 307.1 | 140.1 |
| Paclobutrazol | 3.6 | 294.3 | 70.1 | 125.1 |
| Malathion | 3.6 | 331.0 | 127.2 | 285.2 |
| Myclobutanil | 3.7 | 289.1 | 70.1 | 125.1 |

| Name | Retention Time (min) | Precursor Ion | Product Ion 1 | Product Ion 2 |
|-----------------------------|----------------------|---------------|---------------|---------------|
| Bifenazate | 3.7 | 301.0 | 198.1 | 170.2 |
| Ochratoxin A | 3.8 | 404.2 | 239.1 | 358.3 |
| Fenhexamid | 3.9 | 302.1 | 97.1 | 55.2 |
| Spirotetramat | 4.0 | 374.2 | 302.1 | 216.1 |
| Ethoprophos | 4.1 | 243.1 | 131.1 | 97.1 |
| Fipronil* | 4.1 | 436.8 | 331.8 | 251.9 |
| Fenoxycarb | 4.2 | 302.1 | 88.1 | 116.1 |
| Kresoxim methyl | 4.4 | 314.2 | 267.2 | 222.2 |
| Tebuconazole | 4.4 | 308.1 | 70.1 | 125.1 |
| Diazinon-D10 | 4.6 | 315.2 | 170.2 | - |
| Spinosad (spinosyn A) | 4.6 | 732.4 | 142.2 | 98.1 |
| Diazinon | 4.6 | 305.1 | 169.2 | 153.2 |
| Coumaphos | 4.7 | 363.1 | 227.1 | 307.1 |
| Pyridaben | 4.7 | 365.1 | 309.2 | 147.2 |
| Propiconazole | 4.7 | 342.0 | 159.0 | 69.2 |
| Clofentezine | 4.8 | 303.0 | 138.1 | 102.1 |
| Spinosad (spinosyn D) | 5.0 | 746.5 | 142.3 | 98.4 |
| Spinetoram (spinosyn J) | 5.1 | 748.5 | 142.3 | 98.3 |
| Trifloxystrobin | 5.3 | 409.2 | 186.1 | 145.1 |
| Prallethrin | 5.3 | 301.2 | 123.2 | 105.2 |
| Pyrethrin II | 5.5 | 373.1 | 161.1 | 133.2 |
| Spinetoram (spinosyn L) | 5.6 | 760.5 | 142.2 | 98.1 |
| Piperonyl butoxide | 6.0 | 356.3 | 177.2 | 119.2 |
| Chlorpyrifos | 6.1 | 349.9 | 198.0 | 97.1 |
| Hexythiazox | 6.2 | 353.1 | 228.1 | 168.1 |
| Etoazoxole | 6.6 | 360.2 | 141.1 | 304.2 |
| Spiromesifen | 6.7 | 273.2 | 255.2 | 187.2 |
| Pyrethrin I | 6.9 | 329.2 | 161.2 | 105.2 |
| Cyfluthrin (qualifier) | 6.9 | 453.1 | 193.2 | - |
| Cyfluthrin | 6.9 | 451.1 | 191.2 | - |
| Cypermethrin | 7.1 | 433.1 | 191.0 | 416.0 |
| Fenpyroximate | 7.1 | 422.2 | 366.1 | 138.1 |
| Permethrin- <i>trans</i> | 7.6 | 408.3 | 183.2 | 355.1 |
| Permethrin- <i>cis</i> | 7.9 | 408.3 | 183.2 | 355.1 |
| Avermectin B1a | 7.9 | 890.5 | 305.4 | 567.4 |
| Etofenprox | 8.0 | 394.3 | 177.2 | 359.3 |
| Bifenthrin | 8.2 | 440.0 | 181.2 | 166.2 |
| Acequinocyl precursor ion 1 | 9.4 | 402.3 | 343.2 | 189.0 |
| Acequinocyl precursor ion 2 | 9.4 | 386.0 | 344.2 | 189.1 |

*Analyzed in negative mode.

Table VII: GC-MS/MS transitions.

| Name | Retention Time (min) | Ion Polarity | Precursor Ion | Product Ion |
|--------------------------------|----------------------|--------------|---------------|-------------|
| Atrazine-D5 (IS) (Quan) | 6.82 | Positive | 220.0 | 58.0 |
| Atrazine-D5 (IS) (Qual) | 6.82 | Positive | 205.0 | 127.0 |
| Diazinon-D10 (Quan) | 7.01 | Positive | 183.0 | 139.0 |
| Diazinon-D10 (Qual) | 7.01 | Positive | 183.0 | 168.0 |
| Quintozene (PCNB) (Quan) | 7.03 | Positive | 294.9 | 236.9 |
| Quintozene (PCNB) (Qual) | 7.03 | Positive | 236.8 | 118.9 |
| Methyl parathion (Quan) | 7.50 | Positive | 263.0 | 109.0 |
| Methyl parathion (Qual) | 7.50 | Positive | 263.0 | 79.0 |
| Captan (Quan) | 8.37 | Positive | 184.0 | 149.1 |
| Captan (Qual) | 8.37 | Positive | 184.0 | 134.1 |
| <i>trans</i> -Chlordane (Quan) | 8.41 | Positive | 271.9 | 237.0 |
| <i>trans</i> -Chlordane (Qual) | 8.41 | Positive | 372.9 | 265.9 |
| <i>cis</i> -Chlordane (Quan) | 8.53 | Positive | 372.9 | 265.9 |
| <i>cis</i> -Chlordane (Qual) | 8.53 | Positive | 271.9 | 237.0 |
| Chlorfenapyr (Quan) | 8.80 | Positive | 247.1 | 227.1 |
| Chlorfenapyr (Qual) | 8.80 | Positive | 59.1 | 31.1 |
| Cyfluthrin (Quan) | 10.61 | Positive | 226.0 | 206.0 |
| Cyfluthrin (Qual) | 10.61 | Positive | 163.0 | 127.0 |
| Cypermethrin (Quan) | 10.87 | Positive | 163.0 | 127.1 |
| Cypermethrin (Qual) | 10.87 | Positive | 181.1 | 152.1 |

Table VIII: Cannabinoid Retention Times.

| Compound | Retention Time (min) |
|--|----------------------|
| Cannabidiolic acid (CBDA) | 2.142 |
| Cannabigerol (CBG) | 2.405 |
| Cannabidiol (CBD) | 2.535 |
| Cannabinol (CBN) | 3.776 |
| Delta-9-tetrahydrocannabinol (Delta-9-THC) | 4.753 |
| Tetrahydrocannabinolic acid (THCA) | 6.279 |

Results and Discussion

Method Optimization

The diversity of cannabis matrices and differences in the chemical characteristics of the target analytes requires multiple strategies to ensure accurate results. Analysis of pesticides, mycotoxins, and cannabinoids in cannabis gummies, as required by the state of California, needed a completely different set of experimental conditions from the ones used for brownies in our previous technical article [4]. First, we observed that dealing with pulverized gummies (grinding was conducted using dry ice and a food processor) can be extremely challenging as the matrix becomes very sticky once it gets to room temperature. For that reason, chopping gummies in small pieces was found to be a much better alternative for ease of handling when weighing the desired amount of sample.

The next challenge was finding the best approach to dissolve the matrix prior to extraction in order to obtain reliable data. Since the sample was chopped in small pieces of random size, having a homogeneous matrix was crucial to guarantee satisfactory method reproducibility. First, we evaluated the use of solvents such as acetonitrile and methanol to dissolve the gummy pieces, and to extract analytes of interest as in a simple solvent extraction. However, we found that it was very difficult to solubilize the matrix under those conditions. We also evaluated the use of DMSO for this purpose. Although all gummy pieces were easily dissolved in this solvent, daminozide was not detected in any of the extracts due to ionization suppression caused by the small percentage of DMSO that remained in the injected samples. After these tests, we concluded that hydrating the sample with 5 mL of water followed by vigorous vortexing was the best way to obtain a homogeneous sample. Subsequently, to extract all the contaminants from the matrix, 5 mL of acetonitrile acidified at 1% with acetic acid was added to the dissolved sample.

To separate the organic layer from the aqueous layer, three different Q-sep QuEChERS extraction salts were compared: AOAC (cat.# 25851), unbuffered (cat.# 25847), and EN salts (cat.# 25849). EN salts resulted in the best performance with all the compounds showing recoveries above 83% except for daminozide and its deuterated analogue, which displayed recoveries of around 25%. As for the cleanup step, the effect of four different Q-sep QuEChERS dSPE sorbent mixes (cat.#s 26215, 26216, 26217, 26242) was assessed for all the LC-MS-amenable pesticides. It was confirmed that all the mixes that contained 25 mg of PSA (cat.#s 26215, 26216, 26217) led to significant losses of daminozide and ochratoxin A with only 40% of the original amount present in solution being recovered. The dSPE mix that contained MgSO₄ and C18 (cat.# 26242) did not cause significant losses of pesticides; however, we decided to evaluate the suitability of using the organic extract without any further cleanup step for the quantification of LC-amenable analytes.

To assess the feasibility of analyzing the extracts directly, experiments to investigate absolute matrix effects in gummy extracts without any cleanup step were conducted using the methodology proposed by Matuszowski et al. [5]. Extracts obtained from blank samples were spiked at 5, 15, and 50 ppb final concentrations, and their responses were compared to neat solvent spiked at the same concentration levels using LC-MS/MS. At 5 ppb, 12 pesticides showed matrix effects greater than 120%, and at 15 ppb and 50 ppb, only daminozide showed significant enhancement (Table IX). Based on this and on the poor recovery of daminozide, the use of daminozide-d₆ as the internal standard was crucial to obtaining reliable data.

In regard to the recoveries of GC-amenable pesticides, data corresponding to the evaluation of three dSPE sorbent mixes (cat.#s 26215, 26216, and 26217) showed that, in all the cases, recoveries were above 96% when comparing the response of the cleaned extract vs. the response of the original extract. Considering the high content of sugar and pigments present in gummies extracts, the dSPE sorbent containing PSA, GCB, and magnesium sulfate (cat.# 26217) was chosen for the final sample preparation procedure.

Method Verification

Table X presents results corresponding to limits of quantitation, linearity, accuracy, and precision for the California list of pesticides and mycotoxins determined in gummy matrix. For all the contaminants analyzed via LC-MS/MS, calibration curves were plotted using analyte/internal standard response ratios and a weighing factor of 1/x. For the GC-amenable analytes, only the calibration curve of PCNB was plotted using the analyte/internal standard ratio with diazinon-d₁₀ being chosen as internal standard. Quantification for the rest of the GC-amenable compounds was carried out with external calibration curves (area vs. spiked concentration) because this provided better results than when the internal standard was used. RSDs values below 24% were obtained for all the analytes at all the concentration levels tested. Accuracy values were within 75–118%, and coefficients of determination (R²) were all above 0.99.

Finally, the results for cannabinoids analysis demonstrated that the extract collected for contaminants determination is also suitable for potency testing. Table XI presents data corresponding to the calibration curves prepared in solvent that were used for the quantitation of each cannabinoid. As shown in Table XII, gummy samples spiked with six cannabinoids at 0.2 mg/g exhibited recoveries ranging from 99 to 107%, whereas samples spiked at 0.5 mg/g showed recoveries from 99 to 106%. Representative chromatograms are presented in Figures 2–4.

Table IX: Absolute Matrix Effects (ME) for Pesticides and Mycotoxins in Cannabis Gummies.

| | ME at 5 ppb (%) | RSD | ME at 15 ppb (%) | RSD | ME at 50 ppb (%) | RSD |
|---------------------|--------------------|-----|---------------------|-----|---------------------|-----|
| Daminozide | 216 | 3 | 251 | 7 | 185 | 3 |
| Acephate | 95 | 6 | 85 | 1 | 89 | 6 |
| Oxamyl | 105 | 5 | 96 | 2 | 98 | 5 |
| Flonicamid | 91 | 32 | 90 | 22 | 97 | 17 |
| Methomyl | 104 | 4 | 92 | 3 | 99 | 4 |
| Thiamethoxam | 105 | 5 | 92 | 5 | 96 | 5 |
| Imidacloprid | 114 | 9 | 87 | 9 | 100 | 9 |
| Mevinphos | 102 | 6 | 93 | 3 | 93 | 6 |
| Acetamiprid | 100 | 2 | 88 | 0 | 91 | 2 |
| Dimethoate | 97 | 3 | 92 | 3 | 93 | 3 |
| Thiacloprid | 109 | 7 | 93 | 1 | 93 | 7 |
| Aflatoxin G2 | 102 | 6 | 86 | 0 | 93 | 6 |
| Aflatoxin G1 | 106 | 7 | 93 | 1 | 91 | 7 |
| Aldicarb | 81 | 22 | 85 | 14 | 90 | 22 |
| Aflatoxin B2 | 113 | 11 | 79 | 9 | 95 | 11 |
| Dichlorvos | 126 | 15 | 90 | 6 | 92 | 15 |
| Aflatoxin B1 | 104 | 9 | 93 | 5 | 95 | 9 |
| Imazalil | 92 | 1 | 96 | 3 | 98 | 1 |
| Carbofuran | 106 | 3 | 98 | 1 | 102 | 3 |
| Propoxur | 104 | 3 | 96 | 4 | 96 | 3 |
| Carbaryl | 103 | 4 | 99 | 8 | 93 | 4 |
| Naled | 101 | 4 | 94 | 2 | 84 | 4 |
| Metalaxyl | 105 | 3 | 97 | 2 | 97 | 3 |
| Spiroxamine | 105 | 4 | 94 | 1 | 97 | 4 |
| Chlorantraniliprole | 114 | 6 | 88 | 5 | 100 | 6 |
| Phosmet | 108 | 10 | 100 | 6 | 92 | 10 |
| Azoxystrobin | 104 | 1 | 94 | 1 | 94 | 1 |
| Fludioxonil | 101 | 7 | 90 | 11 | 94 | 7 |
| Methiocarb | 102 | 5 | 94 | 3 | 95 | 5 |
| Dimethomorph | 108 | 15 | 93 | 0 | 93 | 15 |
| Boscalid | 119 | 11 | 94 | 4 | 85 | 11 |
| Paclobutrazol | 110 | 10 | 90 | 2 | 92 | 10 |
| Malathion | 101 | 3 | 92 | 4 | 95 | 3 |
| Myclobutanil | 92 | 11 | 93 | 0 | 95 | 11 |
| Bifenazate | 105 | 2 | 98 | 4 | 99 | 2 |

| | ME at 5 ppb (%) | RSD | ME at 15 ppb (%) | RSD | ME at 50 ppb (%) | RSD |
|--------------------------|--------------------|-----|---------------------|-----|---------------------|-----|
| Ochratoxin A | 111 | 10 | 111 | 0 | 98 | 10 |
| Fenhexamid | 135 | 6 | 102 | 0 | 91 | 6 |
| Spirotetramat | 107 | 7 | 87 | 5 | 95 | 7 |
| Ethoprophos | 104 | 5 | 97 | 2 | 96 | 5 |
| Fipronil | 99 | 12 | 93 | 9 | 96 | 12 |
| Fenoxycarb | 108 | 9 | 97 | 1 | 93 | 9 |
| Kresoxim methyl | 93 | 21 | 108 | 3 | 90 | 21 |
| Tebuconazole | 102 | 6 | 95 | 3 | 97 | 6 |
| Spinosyn A | 100 | 5 | 98 | 2 | 95 | 5 |
| Diazinon | 107 | 2 | 95 | 1 | 96 | 2 |
| Coumaphos | 114 | 7 | 94 | 3 | 97 | 7 |
| Pyridaben | 139 | 18 | 103 | 15 | 94 | 18 |
| Propiconazole | 106 | 2 | 97 | 1 | 95 | 2 |
| Clofentezine | 126 | 14 | 87 | 3 | 92 | 14 |
| Spinosyn D | 104 | 8 | 93 | 8 | 97 | 8 |
| Spinosyn J | 110 | 6 | 92 | 6 | 101 | 6 |
| Trifloxystrobin | 103 | 2 | 95 | 2 | 101 | 2 |
| Prallethrin | 101 | 17 | 103 | 0 | 98 | 17 |
| Pyrethrin II | 92 | 23 | 82 | 5 | 101 | 23 |
| Spinosyn L | 106 | 6 | 95 | 0 | 98 | 6 |
| Piperonyl butoxide | 103 | 2 | 97 | 0 | 96 | 2 |
| Chlorpyrifos | 105 | 10 | 93 | 0 | 91 | 10 |
| Hexythiazox | 121 | 15 | 96 | 0 | 86 | 15 |
| Etoazazole | 102 | 1 | 95 | 1 | 95 | 1 |
| Spiromesifen | 112 | 6 | 95 | 5 | 100 | 6 |
| Pyrethrin I | 131 | 13 | 97 | 4 | 95 | 13 |
| Cyfluthrin | - | - | - | - | 80 | 8 |
| Cypermethrin | - | - | 95 | 16 | 96 | 16 |
| Fenpyroximate | 107 | 6 | 95 | 4 | 94 | 6 |
| Permethrin- <i>trans</i> | 136 | 22 | 103 | 4 | 94 | 22 |
| Permethrin- <i>cis</i> | 113 | 12 | 98 | 8 | 101 | 12 |
| Avermectin B1a | 122 | 17 | 96 | 0 | 87 | 17 |
| Etofenprox | 110 | 4 | 100 | 2 | 101 | 4 |
| Bifenthrin | 135 | 5 | 98 | 6 | 85 | 5 |
| Acequinocyl | 129 | 5 | 95 | 5 | 78 | 5 |

Table X: LOQ, Linearity, Accuracy, and Precision for Pesticides and Mycotoxins in Cannabis Gummies.

| Contaminant | Action level (ng/g) | LOQ (ng/g) | R2 | 10 ng/g (n=4) | | 50 ng/g (n=4) | | 100 ng/g (n=4) | | 500 ng/g (n=4) | |
|---|---------------------|------------|--------|---------------|-----------------|---------------|-----------------|----------------|-----------------|----------------|-----------------|
| | | | | Accuracy (%) | Precision (RSD) | Accuracy (%) | Precision (RSD) | Accuracy (%) | Precision (RSD) | Accuracy (%) | Precision (RSD) |
| Daminozide* | <LOD | 20 | 0.9999 | - | - | 114 | 6 | 114 | 4 | 116 | 4 |
| Acephate | 5000 | 5 | 0.9996 | 100 | 4 | 93 | 3 | 90 | 2 | 88 | 2 |
| Oxamyl | 200 | 5 | 0.999 | 108 | 1 | 102 | 3 | 105 | 4 | 99 | 3 |
| Flonicamid | 2000 | 50 | 0.999 | - | - | 117 | 6 | 105 | 11 | 101 | 4 |
| Methomyl | 100 | 20 | 0.9989 | - | - | 102 | 3 | 103 | 3 | 100 | 1 |
| Thiamethoxam | 4500 | 10 | 0.9988 | 112 | 23 | 108 | 5 | 108 | 6 | 100 | 3 |
| Imidacloprid | 3000 | 10 | 0.9985 | 105 | 17 | 107 | 3 | 109 | 4 | 102 | 5 |
| Mevinphos (I and II)* | <LOD | 20 | 0.9981 | - | - | 101 | 4 | 104 | 7 | 101 | 3 |
| Acetamiprid | 5000 | 10 | 0.9968 | 100 | 8 | 108 | 6 | 109 | 5 | 105 | 1 |
| Dimethoate* | <LOD | 5 | 0.9994 | 109 | 15 | 104 | 3 | 101 | 5 | 99 | 2 |
| Thiacloprid* | <LOD | 20 | 0.9981 | - | - | 100 | 1 | 103 | 4 | 103 | 3 |
| Aflatoxin G2 | 20 [#] | 5 | 0.9957 | 112 | 17 | 101 | 6 | 97 | 2 | - | - |
| Aflatoxin G1 | 20 [#] | 5 | 0.9984 | 114 | 9 | 98 | 1 | 100 | 4 | - | - |
| Aldicarb* | <LOD | 20 | 0.9971 | - | - | 91 | 17 | 104 | 8 | 97 | 4 |
| Aflatoxin B2 | 20 [#] | 5 | 0.9973 | 97 | 23 | 107 | 6 | 94 | 7 | - | - |
| Dichlorvos* | <LOD | 10 | 0.9984 | 98 | 17 | 103 | 4 | 97 | 18 | 106 | 4 |
| Aflatoxin B1 | 20 [#] | 5 | 0.9978 | 113 | 5 | 101 | 6 | 96 | 5 | - | - |
| Imazalil* | <LOD | 5 | 0.9977 | 97 | 19 | 109 | 5 | 107 | 4 | 105 | 3 |
| Carbofuran* | <LOD | 5 | 0.9973 | 93 | 7 | 108 | 1 | 109 | 5 | 99 | 4 |
| Propoxur* | <LOD | 5 | 0.9977 | 108 | 6 | 108 | 2 | 107 | 4 | 102 | 3 |
| Carbaryl | 500 | 5 | 0.9988 | 95 | 14 | 109 | 6 | 108 | 5 | 102 | 2 |
| Naled | 500 | 5 | 0.9968 | 98 | 4 | 112 | 8 | 110 | 3 | 101 | 5 |
| Metalaxyl | 15,000 | 5 | 0.9988 | 101 | 5 | 105 | 4 | 106 | 5 | 99 | 3 |
| Spiroxamine* | <LOD | 5 | 0.9977 | 104 | 6 | 106 | 2 | 105 | 2 | 101 | 3 |
| Chlorantraniliprole | 40,000 | 20 | 0.9971 | - | - | 93 | 7 | 104 | 6 | 104 | 5 |
| Phosmet | 200 | 5 | 0.9992 | 109 | 14 | 105 | 3 | 104 | 3 | 100 | 4 |
| Azoxystrobin | 40,000 | 5 | 0.9992 | 100 | 3 | 104 | 1 | 105 | 4 | 102 | 3 |
| Fludioxonil | 30,000 | 20 | 0.9949 | - | - | 109 | 13 | 97 | 6 | 101 | 6 |
| Methiocarb* | <LOD | 5 | 0.9988 | 116 | 16 | 105 | 5 | 107 | 5 | 101 | 4 |
| Dimethomorph (I and II) | 20,000 | 10 | 0.999 | 75 | 14 | 101 | 6 | 93 | 8 | 101 | 5 |
| Boscalid | 10,000 | 10 | 0.9964 | 108 | 15 | 102 | 12 | 102 | 2 | 103 | 3 |
| Paclobutrazol* | <LOD | 10 | 0.9979 | 99 | 9 | 106 | 1 | 108 | 3 | 100 | 4 |
| Malathion | 5000 | 10 | 0.9989 | 117 | 12 | 112 | 3 | 106 | 2 | 101 | 3 |
| Myclobutanil | 9000 | 10 | 0.9986 | 102 | 22 | 102 | 6 | 104 | 3 | 101 | 3 |
| Bifenazate | 5000 | 10 | 0.9994 | 118 | 20 | 110 | 3 | 104 | 10 | 102 | 5 |
| Ochratoxin A | 20 | 10 | 0.9957 | 99 | 8 | 99 | 24 | 104 | 8 | - | - |
| Fenhexamid | 10,000 | 10 | 0.9969 | 96 | 21 | 109 | 4 | 106 | 3 | 107 | 5 |
| Spirotetramat | 13,000 | 10 | 0.9987 | 83 | 18 | 107 | 6 | 106 | 3 | 105 | 3 |
| Ethoprophos* | <LOD | 5 | 0.9985 | 101 | 3 | 106 | 4 | 104 | 1 | 102 | 3 |
| Fipronil* | <LOD | 20 | 0.998 | - | - | 96 | 5 | 104 | 8 | 103 | 4 |
| Fenoxycarb* | <LOD | 10 | 0.9967 | 115 | 5 | 106 | 4 | 107 | 2 | 102 | 3 |
| Kresoxym-methyl | 1000 | 10 | 0.9993 | 112 | 20 | 104 | 10 | 102 | 4 | 103 | 3 |
| Tebuconazole | 2000 | 5 | 0.999 | 100 | 3 | 110 | 1 | 105 | 3 | 101 | 4 |
| Spinosad - spinosyn A (71 %) ^a | 3000 ^y | 7.1 | 0.9988 | 117 | 2 | 110 | 3 | 108 | 2 | 102 | 1 |
| Diazinon | 200 | 5 | 0.9997 | 104 | 1 | 102 | 1 | 104 | 2 | 101 | 1 |

Continued

Table X (cont.)

| Contaminant | Action level (ng/g) | LOQ (ng/g) | R2 | 10 ng/g (n=4) | | 50 ng/g (n=4) | | 100 ng/g (n=4) | | 500 ng/g (n=4) | |
|--|---------------------|------------|--------|---------------|-----------------|---------------|-----------------|----------------|-----------------|----------------|-----------------|
| | | | | Accuracy (%) | Precision (RSD) | Accuracy (%) | Precision (RSD) | Accuracy (%) | Precision (RSD) | Accuracy (%) | Precision (RSD) |
| Coumaphos* | <LOD | 10 | 0.9991 | 107 | 11 | 109 | 5 | 108 | 2 | 103 | 2 |
| Pyridaben | 3000 | 50 | 0.9994 | - | - | 93 | 12 | 113 | 6 | 103 | 1 |
| Propiconazole | 20,000 | 5 | 0.9991 | 96 | 6 | 108 | 4 | 104 | 5 | 103 | 3 |
| Clofentezine | 500 | 20 | 0.9978 | - | - | 102 | 5 | 105 | 3 | 106 | 4 |
| Spinosad - spinosyn D (29%) ^b | 3000 ^g | 2.9 | 0.9993 | 106 | 8 | 101 | 2 | 103 | 5 | 104 | 4 |
| Spinetoram - spinosyn J (80%) ^c | 3000 ^g | 4 | 0.9995 | 104 | 7 | 102 | 3 | 108 | 4 | 102 | 2 |
| Trifloxystrobin | 30,000 | 5 | 0.9996 | 107 | 3 | 106 | 3 | 106 | 3 | 104 | 2 |
| Prallethrin | 400 | 10 | 0.9967 | 107 | 21 | 99 | 16 | 111 | 6 | 101 | 3 |
| Pyrethrin II (34%) ^f | 1000 ^g | 17 | 0.9977 | - | - | 94 | 11 | 112 | 4 | 106 | 6 |
| Spinetoram - spinosyn L (20%) ^d | 3000 ^g | 2 | 0.9993 | 107 | 10 | 109 | 4 | 107 | 3 | 102 | 1 |
| Piperonyl Butoxide | 8000 | 5 | 0.9998 | 110 | 5 | 99 | 2 | 94 | 3 | 95 | 5 |
| Chlorpyrifos* | <LOD | 20 | 0.9995 | - | - | 103 | 9 | 106 | 4 | 104 | 5 |
| Hexythiazox | 2000 | 10 | 0.9975 | 104 | 15 | 102 | 3 | 107 | 4 | 106 | 5 |
| Etoazole | 1500 | 5 | 0.9995 | 103 | 4 | 104 | 2 | 103 | 1 | 100 | 1 |
| Spiromesifen | 12,000 | 5 | 0.9992 | 98 | 8 | 109 | 7 | 111 | 2 | 102 | 2 |
| Pyrethrin I (54%) ^e | 1000 ^g | 11 | 0.9977 | - | - | 98 | 5 | 105 | 12 | 104 | 5 |
| Cyfluthrin | 1000 | 50 | 0.999 | - | - | 93 | 11 | 102 | 23 | 115 | 10 |
| Cypermethrin | 1000 | 50 | 0.9961 | - | - | 115 | 11 | 98 | 18 | 104 | 6 |
| Fenpyroximate | 2000 | 5 | 0.9995 | 102 | 9 | 109 | 3 | 109 | 2 | 105 | 2 |
| Permethrin-trans (59%) ^h | 20,000 ^c | 12 | 0.9994 | - | - | 99 | 8 | 104 | 5 | 101 | 3 |
| Permethrin-cis (41%) ^g | 20,000 ^c | 8 | 0.9996 | 95 | 5 | 104 | 5 | 103 | 4 | 101 | 3 |
| Avermectin B1a | 300 | 50 | 0.9988 | - | - | 114 | 3 | 108 | 2 | 105 | 2 |
| Etofenprox* | <LOD | 5 | 0.9994 | 104 | 7 | 107 | 1 | 106 | 2 | 103 | 1 |
| Bifenthrin | 500 | 5 | 0.999 | 99 | 5 | 103 | 2 | 108 | 6 | 104 | 2 |
| Acequinocyl | 4000 | 10 | 0.9997 | 104 | 7 | 109 | 3 | 108 | 3 | 106 | 2 |
| Quintozene (PCNB) (GC) | 200 | 10 | 0.9966 | 110 | 14 | 102 | 5 | 101 | 1 | 97 | 3 |
| Methyl parathion (GC)* | <LOD | 5 | 0.9934 | 89 | 10 | 89 | 5 | 87 | 3 | 89 | 6 |
| Captan (GC) | 5000 | 10 | 0.9924 | 110 | 14 | 96 | 10 | 94 | 18 | 92 | 7 |
| Chlordane (GC)* | <LOD | 20 | 0.9913 | - | - | 105 | 9 | 93 | 8 | 85 | 10 |
| Chlorfenapyr (GC)* | <LOD | 10 | 0.9924 | 97 | 9 | 90 | 19 | 89 | 6 | 85 | 13 |
| Cyfluthrin (GC) | 1000 | 5 | 0.9935 | 107 | 10 | 92 | 18 | 91 | 7 | 89 | 11 |
| Cypermethrin (GC) | 1000 | 5 | 0.9938 | 95 | 9 | 83 | 17 | 97 | 15 | 93 | 9 |

*Category I pesticides, LOQ ≤100 ng/g

^a Spinosad - spinosyn A: Conc. 1: 7 ng/g; Conc. 2: 35.5 ng/g; Conc. 3: 71 ng/g; Conc. 4: 355 ng/g

^b Spinosad - spinosyn D: Conc. 1: 3 ng/g; Conc. 2: 14.5 ng/g; Conc. 3: 29 ng/g; Conc. 4: 145 ng/g

^c Spinetoram - spinosyn J: Conc. 1: 8 ng/g; Conc. 2: 40 ng/g; Conc. 3: 80 ng/g; Conc. 4: 400 ng/g

^d Spinetoram - spinosyn L: Conc. 1: 2 ng/g; Conc. 2: 10 ng/g; Conc. 3: 20 ng/g; Conc. 4: 100 ng/g

^e Pyrethrin I: Conc. 1: 5 ng/g; Conc. 2: 27 ng/g; Conc. 3: 54 ng/g; Conc. 4: 270 ng/g

^f Pyrethrin II: Conc. 1: 3 ng/g; Conc. 2: 17 ng/g; Conc. 3: 34 ng/g; Conc. 4: 170 ng/g

^g Permethrin-cis: Conc. 1: 4 ng/g; Conc. 2: 20.5 ng/g; Conc. 3: 41 ng/g; Conc. 4: 205 ng/g

^h Permethrin-trans: Conc. 1: 6 ng/g; Conc. 2: 29.5 ng/g; Conc. 3: 59 ng/g; Conc. 4: 295 ng/g

ⁱ Total of aflatoxin B1, B2, G1, and G2 should not exceed 20 ng/g.

^j Total spinosad should not exceed 3000 ng/g.

^k Total spinetoram should not exceed 3000 ng/g.

^l Total pyrethrins should not exceed 1000 ng/g.

^m Total permethrins should not exceed 20,000 ng/g.

Table XI: Linearity for Cannabinoids in Cannabis Gummies.


| Cannabinoids | Retention time | R ² | Equation |
|--|----------------|----------------|------------------------------|
| Cannabidiolic acid (CBDA) | 2.142 | 0.9993 | $y = 2.08e+004x + 1.03e+003$ |
| Cannabigerol (CBG) | 2.405 | 0.9981 | $y = 1.16e+004x + 1.85e+003$ |
| Cannabidiol (CBD) | 2.535 | 0.9972 | $y = 1.17e+004x + 1.42e+003$ |
| Cannabinol (CBN) | 3.776 | 0.9980 | $y = 2.70e+004x + 6.14e+004$ |
| Delta-9-tetrahydrocannabinol (Delta-9-THC) | 4.753 | 0.9970 | $y = 1.06e+004x + 7.29e+003$ |
| Tetrahydrocannabinolic acid (THCA) | 6.279 | 0.9986 | $y = 1.78e+004x - 1.28e+003$ |

Table XII: Accuracy and Precision for Cannabinoids in Cannabis Gummies.

| Cannabinoid/Spike Level | Diluted Extract Concentration (ppm) | | Average (ppm) | SD | RSD (%) | Undiluted Extract Conc. (ppm) | Estimated Sample Concentration (mg/g) | Accuracy (%) | Percent Error |
|--|-------------------------------------|-------------|---------------|-----|---------|-------------------------------|---------------------------------------|--------------|---------------|
| Gummy spiked at 0.2 mg/g | Replicate 1 | Replicate 2 | | | | | | | |
| Cannabidiolic acid (CBDA) | 4.2 | 4.4 | 4.3 | 0.1 | 3 | 43 | 0.2 | 107 | 7 |
| Cannabigerol (CBG) | 4.0 | 3.9 | 4.0 | 0.1 | 2 | 40 | 0.2 | 99 | 1 |
| Cannabidiol (CBD) | 4.1 | 4.2 | 4.1 | 0.1 | 2 | 41 | 0.2 | 103 | 3 |
| Cannabinol (CBN) | 4.0 | 4.1 | 4.1 | 0.1 | 3 | 41 | 0.2 | 101 | 1 |
| Delta-9 tetrahydrocannabinol (Delta 9 THC) | 4.0 | 4.2 | 4.1 | 0.1 | 3 | 41 | 0.2 | 103 | 3 |
| Tetrahydrocannabinolic acid (THCA) | 4.1 | 4.3 | 4.2 | 0.1 | 3 | 42 | 0.2 | 105 | 5 |
| | | | | | | | | | |
| Gummy spiked at 0.5 mg/g | | | | | | | | | |
| Cannabidiolic acid (CBDA) | 10.4 | 10.7 | 10.6 | 0.2 | 2 | 106 | 0.5 | 106 | 6 |
| Cannabigerol (CBG) | 9.5 | 10.3 | 9.9 | 0.5 | 5 | 99 | 0.5 | 99 | 1 |
| Cannabidiol (CBD) | 9.8 | 10.6 | 10.2 | 0.5 | 5 | 102 | 0.5 | 102 | 2 |
| Cannabinol (CBN) | 9.8 | 10.3 | 10.0 | 0.4 | 4 | 100 | 0.5 | 100 | 0 |
| Delta-9-tetrahydrocannabinol (Delta-9-THC) | 9.9 | 10.3 | 10.1 | 0.3 | 3 | 101 | 0.5 | 101 | 1 |
| Tetrahydrocannabinolic acid (THCA) | 10.3 | 9.9 | 10.1 | 0.2 | 2 | 101 | 0.5 | 101 | 1 |

Figure 2: LC-MS/MS chromatogram of an extract obtained from blank gummy spiked with pesticides and mycotoxins at 100 ng/g.


- LC amenable pesticides and mycotoxins regulated by California (2020).
- Quantify low ng/g concentrations in cannabis edibles.



| Peaks | t_r (min) | Precursor Ion | Product Ion 1 | Product Ion 2 | Polarity | Peaks | t_r (min) | Precursor Ion | Product Ion 1 | Product Ion 2 | Polarity |
|-----------------------------|----------------|------------------|------------------|------------------|----------|---|----------------|------------------|------------------|------------------|----------|
| 1. Daminozide-d6 | 0.7 | 167.0 | 149.3 | 49.3 | + | 63. Spinosyn L (Spinetoram) | 5.6 | 760.5 | 142.2 | 98.1 | + |
| 2. Daminozide | 0.7 | 161.1 | 44.1 | 143.2 | + | 64. Piperonyl butoxide | 6.0 | 356.3 | 177.2 | 119.2 | + |
| 3. Acephate | 1.7 | 184.0 | 143.1 | 95.1 | + | 65. Chlorpyrifos | 6.1 | 349.9 | 198.0 | 97.1 | + |
| 4. Oxamyl | 2.0 | 237.1 | 72.1 | 90.1 | + | 66. Hexythiazox | 6.2 | 353.1 | 228.1 | 168.1 | + |
| 5. Flonicamid | 2.1 | 230.1 | 203.1 | 174.1 | + | 67. Etoxazole | 6.6 | 360.2 | 141.1 | 304.2 | + |
| 6. Methomyl | 2.1 | 163.1 | 88.1 | 106.1 | + | 68. Spiromesifen | 6.7 | 273.2 | 255.2 | 187.2 | + |
| 7. Thiamethoxam | 2.1 | 292.0 | 211.1 | 181.1 | + | 69. Pyrethrin I | 6.9 | 329.2 | 161.2 | 105.2 | + |
| 8. Imidacloprid | 2.3 | 256.1 | 209.1 | 175.1 | + | 70. Cyfluthrin (qualifier) | 6.9 | 453.1 | 193.2 | - | + |
| 9. Mevinphos | 2.4 | 225.1 | 127.1 | 193.2 | + | 71. Cyfluthrin | 6.9 | 451.1 | 191.2 | - | + |
| 10. Acetamiprid | 2.4 | 223.0 | 126.1 | 56.1 | + | 72. Cypermethrin | 7.1 | 433.1 | 191.0 | 416.0 | + |
| 11. Dimethoate-d6 | 2.4 | 236.1 | 205.1 | - | + | 73. Fenpyroximate | 7.1 | 422.2 | 366.1 | 138.1 | + |
| 12. Dimethoate | 2.4 | 230.0 | 199.1 | 125.1 | + | 74. trans-Permethrin | 7.6 | 408.3 | 183.2 | 355.1 | + |
| 13. Thiacloprid | 2.5 | 253.0 | 126.0 | 90.1 | + | 75. cis-Permethrin | 7.9 | 408.3 | 183.2 | 355.1 | + |
| 14. Aflatoxin G2 | 2.5 | 331.2 | 189.3 | 115.2 | + | 76. Avermectin B1a | 7.9 | 890.5 | 305.4 | 567.4 | + |
| 15. Aflatoxin G1 | 2.5 | 329.2 | 243.2 | 215.3 | + | 77. Etofenprox | 8.0 | 394.3 | 177.2 | 359.3 | + |
| 16. Aldicarb | 2.6 | 116.0 | 89.2 | 70.2 | + | 78. Bifenthrin | 8.2 | 440.0 | 181.2 | 166.2 | + |
| 17. Aflatoxin B2 | 2.6 | 315.3 | 287.2 | 243.3 | + | 79. Acequinocyl (precursor ion 1) | 9.4 | 402.3 | 343.2 | 189.0 | + |
| 18. Dichlorvos | 2.7 | 220.9 | 109.1 | 79.2 | + | 80. Acequinocyl (precursor ion 2) | 9.4 | 386.0 | 344.2 | 189.1 | + |
| 19. Dichlorvos-d6 | 2.7 | 227.0 | 115.1 | - | + | Column Raptor ARC-18 (cat.# 9314A12) Dimensions: 100 mm x 2.1 mm ID Particle Size: 2.7 µm Pore Size: 90 Å Guard Column: Raptor ARC-18 EXP guard column cartridge 5 mm, 2.1 mm ID, 2.7 µm (cat.# 9314A0252) Temp.: 40 °C Sample California pesticide standard #1 (cat.# 34124); California pesticide standard #2 (cat.# 34125); California pesticide standard #3 (cat.# 34126); California pesticide standard #4 (cat.# 34127); California pesticide standard #5 (cat.# 34128); California pesticide standard #6 (cat.# 34129); Dimethoate-d6 (cat.# 31988); Dichlorvos-d6 (cat.# 31987); Carbaryl-d7 (cat.# 31985); Diazinon-d10 (cat.# 31986); Atrazine-d5 (cat.# 31984); Diuron-d6 (cat.# 31989); Liuron-d6 (cat.# 31990); Aflatoxins standard (cat.# 34121); Ochratoxin A (cat.# 34122); Compounds not present in these mixes were obtained separately. Diluent: 75:25 Acetonitrile:water Conc.: 3.75-15 ng/mL (Expected concentration range in extract of gummy initially spiked at 100 ng/g.) Inj. Vol.: 2 µL Mobile Phase A: Water, 2 mM ammonium formate, 0.1% formic acid B: Methanol, 2 mM ammonium formate, 0.1% formic acid | | | | | |
| 20. Aflatoxin B1 | 2.7 | 313.2 | 241.2 | 128.2 | + | | | | | | |
| 21. Imazalil | 2.7 | 297.0 | 159.0 | 201.0 | + | | | | | | |
| 22. Carbofuran | 2.7 | 222.1 | 123.1 | 165.2 | + | | | | | | |
| 23. Propoxur | 2.7 | 210.1 | 111.1 | 93.1 | + | | | | | | |
| 24. Carbaryl-d7 | 2.8 | 209.2 | 152.2 | - | + | | | | | | |
| 25. Carbaryl | 2.8 | 202.1 | 145.1 | 127.1 | + | | | | | | |
| 26. Diuron-d6 | 3.0 | 239.1 | 78.2 | - | + | | | | | | |
| 27. Atrazine-d5 | 3.0 | 221.2 | 179.1 | - | + | | | | | | |
| 28. Naled | 3.1 | 397.8 | 127.1 | 109.1 | + | | | | | | |
| 29. Metalaxyl | 3.1 | 280.2 | 220.2 | 192.2 | + | | | | | | |
| 30. Spiroxamine | 3.1 | 298.3 | 144.2 | 100.2 | + | | | | | | |
| 31. Chlorantraniliprole | 3.2 | 483.9 | 452.9 | 285.9 | + | | | | | | |
| 32. Phosmet | 3.2 | 318.0 | 160.1 | 77.2 | + | | | | | | |
| 33. Azoxystrobin | 3.3 | 404.0 | 372.1 | 344.1 | + | | | | | | |
| 34. Linuron-d6 | 3.3 | 255.1 | 160.1 | - | + | | | | | | |
| 35. Fludioxonil | 3.4 | 247.0 | 180.0 | 126.0 | - | | | | | | |
| 36. Methiocarb | 3.4 | 226.1 | 169.1 | 121.1 | + | | | | | | |
| 37. Dimethomorph | 3.5 | 388.2 | 301.2 | 165.3 | + | | | | | | |
| 38. Boscalid | 3.5 | 342.9 | 307.1 | 140.1 | + | | | | | | |
| 39. Paclobutrazol | 3.6 | 294.3 | 70.1 | 125.1 | + | | | | | | |
| 40. Malathion | 3.6 | 331.0 | 127.2 | 285.2 | + | | | | | | |
| 41. Myclobutanil | 3.7 | 289.1 | 70.1 | 125.1 | + | | | | | | |
| 42. Bifenazate | 3.7 | 301.0 | 198.1 | 170.2 | + | | | | | | |
| 43. Ochratoxin A | 3.8 | 404.2 | 239.1 | 358.3 | + | | | | | | |
| 44. Fenhexamid | 3.9 | 302.1 | 97.1 | 55.2 | + | | | | | | |
| 45. Spirotetramat | 4.0 | 374.2 | 302.1 | 216.1 | + | | | | | | |
| 46. Ethoprophos | 4.1 | 243.1 | 131.1 | 97.1 | + | | | | | | |
| 47. Fipronil | 4.1 | 436.8 | 331.8 | 251.9 | - | | | | | | |
| 48. Fenoxycarb | 4.2 | 302.1 | 88.1 | 116.1 | + | | | | | | |
| 49. Kresoxim-methyl | 4.4 | 314.2 | 267.2 | 222.2 | + | | | | | | |
| 50. Tebuconazole | 4.4 | 308.1 | 70.1 | 125.1 | + | | | | | | |
| 51. Diazinon-d10 | 4.6 | 315.2 | 170.2 | - | + | | | | | | |
| 52. Spinosyn A (Spinosad) | 4.6 | 732.4 | 142.2 | 98.1 | + | | | | | | |
| 53. Diazinon | 4.6 | 305.1 | 169.2 | 153.2 | + | | | | | | |
| 54. Coumaphos | 4.7 | 363.1 | 227.1 | 307.1 | + | | | | | | |
| 55. Pyridaben | 4.7 | 365.1 | 309.2 | 147.2 | + | | | | | | |
| 56. Propiconazole | 4.7 | 342.0 | 159.0 | 69.2 | + | | | | | | |
| 57. Clofentazine | 4.8 | 303.0 | 138.1 | 102.1 | + | | | | | | |
| 58. Spinosyn D (Spinosad) | 5.0 | 746.5 | 142.3 | 98.4 | + | | | | | | |
| 59. Spinosyn J (Spinetoram) | 5.1 | 748.5 | 142.3 | 98.3 | + | | | | | | |
| 60. Trifloxystrobin | 5.3 | 409.2 | 186.1 | 145.1 | + | | | | | | |
| 61. Prallethrin | 5.3 | 301.2 | 123.2 | 105.2 | + | | | | | | |
| 62. Pyrethrin II | 5.5 | 373.1 | 161.1 | 133.2 | + | | | | | | |
| | | | | | | Detector MS/MS Ion Mode: ESI+/ESI- Mode: MRM Instrument UHPLC Notes Gummies were manually chopped into small pieces, and 1 g of sample was weighed in a 50 mL polypropylene tube. The sample was mixed with 5 mL of water and then vigorously vortexed until all gummy pieces were fully solubilized. The sample was fortified with pesticides and mycotoxins at 100 ng/g. A mix of internal standards was added at 200 ng/g. The spiked sample was further vortexed for 30 sec. 5 mL of acetonitrile acidified with 1% acetic acid was added to the sample, and this was followed by a 30 sec vortex agitation. Then, a pouch of European EN 15662 QuEChERS extraction salts (cat.# 25849) was added to the sample. The sample was vortexed for 30 sec and then centrifuged for 5 min. 750 µL of organic extract was mixed with 250 µL of water. 2 µL of final extract was injected into the LC-MS/MS system. | | | | | |


Figure 3: GC-MS/MS chromatogram of an extract obtained from blank gummy spiked with pesticides and mycotoxins at 100 ng/g.

- GC amenable pesticides regulated by California (2020).
- Quantify pesticides in cannabis edibles at low ng/g concentrations.



| Peaks | t_r (min) | Polarity | Precursor Ion | Product Ion | Transition Type | Column Sample |
|----------------------------------|----------------|----------|------------------|----------------|--------------------|--|
| 1. Atrazine-d5 | 6.82 | Positive | 220.0 | 58.0 | Quantifier | Rxi-5ms, 30 m, 0.25 mm ID, 0.25 μ m (cat.# 13423) California pesticide standard #1 (cat.# 34124) California pesticide standard #2 (cat.# 34125) California pesticide standard #3 (cat.# 34126) California pesticide standard #4 (cat.# 34127) California pesticide standard #5 (cat.# 34128) California pesticide standard #6 (cat.# 34129) Atrazine-d5 (cat.# 31984) Diazinon-d10 (cat.# 31986) Diluent: Acetonitrile Conc.: 2.5-10 ng/mL Expected concentration range in extract after extracting from gummy fortified at 100 ng/g (final extract was diluted in half with acetonitrile). Injection Inj. Vol.: 1 μ L splitless Liner: Topaz 4.0 mm ID single taper inlet liner w/wool (cat.# 23447) Inj. Temp.: 250 °C Purge Flow: 5 mL/min Oven Oven Temp.: 90 °C (hold 1 min) to 310 °C at 25 °C/min (hold 10 min) Carrier Gas He, constant flow Flow Rate: 1.4 mL/min Detector MS/MS Transfer Line Temp.: 290 °C Analyzer Type: Quadrupole Source Temp.: 330 °C Electron Energy: 70 eV Tune Type: PFTBA Ionization Mode: EI Instrument Notes Thermo Scientific TSQ 8000 Triple Quadrupole GC-MS Gummies were manually chopped into small pieces, and 1 g of sample was weighed in a 50 mL polypropylene tube. The sample was mixed with 5 mL of water and then vigor- ously vortexed until all gummy pieces were fully solubilized. The sample was fortified with pesticides and mycotoxins at 100 ng/g. A mix of internal standards was added at 200 ng/g. The spiked sample was further vortexed for 30 sec. 5 mL of acetonitrile acidified with 1% acetic acid was added to the sample, and this was followed by 30 sec vortex agitation. Then, a pouch of European EN 15662 QuEChERS extraction salts (cat.# 25849) was added to the sample. The sample was vortexed for 30 sec and then centrifuged for 5 min. 1.9 mL of supernatant was transferred to a Q-sep QuEChERS dSPE tube containing pre-weighed magnesium sulfate, PSA, and GCB (cat.# 26217). After vortexing and centrifuging, 500 μ L of extract was mixed with 500 μ L of acidified acetonitrile. 1 μ L of final extract was injected into the GC-MS/MS system. |
| 2. Atrazine-d5 | 6.82 | Positive | 205.0 | 127.0 | Qualifier | |
| 3. Diazinon-d10 (diethyl-d10) | 7.01 | Positive | 183.0 | 139.0 | Quantifier | |
| 4. Diazinon-d10 (diethyl-d10) | 7.01 | Positive | 183.0 | 168.0 | Qualifier | |
| 5. Quintozene | 7.03 | Positive | 294.9 | 236.9 | Quantifier | |
| 6. Quintozene | 7.03 | Positive | 236.8 | 118.9 | Qualifier | |
| 7. Methyl parathion | 7.50 | Positive | 263.0 | 109.0 | Quantifier | |
| 8. Methyl parathion | 7.50 | Positive | 263.0 | 79.0 | Qualifier | |
| 9. Captan | 8.37 | Positive | 184.0 | 149.1 | Quantifier | |
| 10. Captan | 8.37 | Positive | 184.0 | 134.1 | Qualifier | |
| 11. trans-Chlordane | 8.41 | Positive | 271.9 | 237.0 | Quantifier | |
| 12. trans-Chlordane | 8.41 | Positive | 372.9 | 265.9 | Qualifier | |
| 13. cis-Chlordane | 8.53 | Positive | 372.9 | 265.9 | Quantifier | |
| 14. cis-Chlordane | 8.53 | Positive | 271.9 | 237.0 | Qualifier | |
| 15. Chlorfenapyr | 8.80 | Positive | 247.1 | 227.1 | Quantifier | |
| 16. Chlorfenapyr | 8.80 | Positive | 59.1 | 31.1 | Qualifier | |
| 17. Cyfluthrin | 10.61 | Positive | 226.0 | 206.0 | Quantifier | |
| 18. Cyfluthrin | 10.61 | Positive | 163.0 | 127.0 | Qualifier | |
| 19. Cypermethrin | 10.87 | Positive | 163.0 | 127.1 | Qualifier | |
| 20. Cypermethrin | 10.87 | Positive | 181.1 | 152.1 | Qualifier | |

Figure 4: HPLC-UV chromatogram of an extract obtained from blank gummy spiked with six cannabinoids at 0.2 mg/g.



Column Raptor ARC-18 (cat.# 9314A62)
Dimensions: 150 mm x 2.1 mm ID
Particle Size: 2.7 μ m
Pore Size: 90 Å
Guard Column: Raptor ARC-18 EXP guard column cartridge 5 mm, 2.1 mm ID, 2.7 μ m (cat.# 9314A0252)
Temp.: 30 °C
Sample Cannabinoids standard (cat.# 34014)
 Cannabigerol (cat.# 34091)
 Δ^9 -Tetrahydrocannabinol (cat.# 34067)
 Δ^9 -Tetrahydrocannabinolic acid A (cat.# 34111)
Diluent: 75:25 Acetonitrile:water
Conc.: Expected concentration of 4 ppm in final extract from gummy initially spiked at 0.2 mg/g.
Inj. Vol.: 2 μ L

Mobile Phase
A: Water, 5 mM ammonium formate, 0.1% formic acid
B: Acetonitrile, 0.1% formic acid

| Time (min) | Flow (mL/min) | %A | %B |
|------------|---------------|----|----|
| 0.00 | 0.4 | 25 | 75 |
| 10.00 | 0.4 | 25 | 75 |

Detector UV/Vis @ 228 nm
Instrument UHPLC
Notes

Gummies were manually chopped into small pieces, and 1 g of sample was weighed in a 50 mL polypropylene tube. The sample was mixed with 5 mL of water and then vigorously vortexed until all gummy pieces were fully solubilized. The sample was fortified with cannabinoids at 0.2 mg/g. The spiked sample was further vortexed for 30 sec. 5 mL of acetonitrile acidified with 1% acetic acid was added to the sample, and this was followed by 30 sec vortex agitation. Then, a pouch of European EN 15662 QuEChERS extraction salts (cat.# 25849) was added to the sample. The sample was vortexed for 30 sec and then centrifuged for 5 min. 100 μ L of organic extract was mixed with 900 μ L of 75:25 acetonitrile:water. 2 μ L of final extract was injected into the HPLC-UV system.

Conclusion

An easy and effective workflow for the analysis of pesticides, mycotoxins, and cannabinoids in cannabis gummies was developed. Sample preparation conditions involved matrix homogenization of gummy pieces with water, extraction of analytes using acidified acetonitrile followed by a salting-out step using Q-sep QuEChERS extraction salts; extract dilution (for LC-MS/MS amenable contaminants and cannabinoids); and dSPE cleanup using magnesium sulfate, PSA, and GCB (for GC-MS/MS amenable pesticides). Satisfactory results in terms of LOQ, linearity, accuracy, and precision were obtained for all the target contaminants. In addition, our data demonstrated that the proposed methodology is suitable for potency testing with accuracy values ranging from 99 to 107% for the six cannabinoids listed in the cannabis regulations of the state of California. Overall, the presented workflow streamlines work for cannabis testing labs by enabling the satisfactory quantitation of multiple analyte classes in gummy samples using a single extract.

References

- [1] Text of Regulations, Bureau of Cannabis Control, California Code of Regulations, <https://cannabis.ca.gov/wp-content/uploads/sites/13/2019/01/Order-of-Adoption-Clean-Version-of-Text.pdf>, (accessed 8 November 2019).
- [2] X. Wang, D. Mackowsky, J. Searfoss, M. J. Telepchak, Determination of cannabinoid content and pesticide residues in cannabis edibles and beverages, Cannabis Sci. and Tech. (2018). <https://www.cannabissciencetech.com/view/determination-cannabinoid-content-and-pesticide-residues-cannabis-edibles-and-beverages>.
- [3] Fast, low-solvent analysis of cannabinoids increases lab productivity and decreases solvent costs, Restek Corporation. https://www.restek.com/Technical-Resources/Technical-Library/Foods-Flavors-Fragrances/fff_FFFA3123-UNV (accessed 22 January 2021).
- [4] N. Reyes-Garcés, C. Myers, Analysis of pesticides and mycotoxins cannabis brownies, Restek Corporation. https://www.restek.com/Technical-Resources/Technical-Library/Foods-Flavors-Fragrances/fff_FFAN3149-UNV (accessed 3 December 2020).
- [5] B. K. Matuszewski, M. L. Constanzer, C. M. Chavez-Eng, Strategies for the assessment of matrix effect in quantitative bioanalytical methods based on HPLC-MS/MS, Anal. Chem. 75 (2003) 3019–3030. <https://pubs.acs.org/doi/10.1021/ac020361s>

Raptor ARC-18 LC Columns (USP L1)

- Ideal for high-throughput LC-MS/MS applications with minimal sample preparation.
- Well-balanced retention profile for better detection and integration of large, multiclass analyte lists.
- Sterically protected to endure low-pH mobile phases without sacrificing retention or peak quality.
- Part of Restek's Raptor LC column line featuring 1.8, 2.7, and 5 μm SPP core-shell silica.



| ID | Length | qty. | cat.# |
|---|--------|------|---------|
| 2.7 μm Particles | | | |
| 2.1 mm | 100 mm | ea. | 9314A12 |
| | 150 mm | ea. | 9314A62 |

Stationary Phase Category: C18, octadecylsilane (L1)

Ligand Type: Sterically protected C18

Particle: 1.8 μm , 2.7 μm , or 5 μm superficially porous silica (SPP or "core-shell")

Pore Size: 90 Å

Carbon Load: 7% (1.8 μm), 7% (2.7 μm), 5% (5 μm)

End-Cap: no

Surface Area: 125 m^2/g (1.8 μm), 130 m^2/g (2.7 μm), or 100 m^2/g (5 μm)

Recommended Usage:

pH Range: 1.0–8.0

Maximum Temperature: 80 °C

Maximum Pressure: 1034 bar/15000 psi* (1.8 μm), 600 bar/8700 psi (2.7 μm); 400 bar/5800 psi (5 μm)

* For maximum lifetime, recommended maximum pressure for 1.8 μm particles is 830 bar/12000 psi.

Properties:

- Well-balanced retention profile.
- Sterically protected and acid resistant to resist harsh, low-pH mobile phases.
- Ideal for use with sensitive detectors like mass spec.

Switch to an ARC-18 when:

- You are analyzing large, multiclass lists by LC-MS/MS.
- Strongly acidic (pH 1–3) mobile phases are required.

Raptor EXP Guard Column Cartridges

- Free-Turn architecture lets you change cartridges by hand without breaking inlet/outlet fluid connections—no tools needed.
- Patented titanium hybrid ferrules can be installed repeatedly without compromising high-pressure seal.
- Auto-adjusting design provides ZDV (zero dead volume) connection to any 10-32 female port.
- Guard column cartridges require EXP direct connect holder (cat.# 25808).
- Pair with EXP hand-tight fitting (cat.# 25937–25938) for tool-free installation.



| Description | Particle Size | Size | qty. | cat.# |
|--|-------------------|------------|-------|-----------|
| Raptor ARC-18 EXP Guard Column Cartridge | 2.7 μm | 5 x 2.1 mm | 3-pk. | 9314A0252 |

Maximum cartridge pressure: 1034 bar/15,000 psi* (UHPLC), 600 bar/8700 psi (2.7 μm); 400 bar/5800 psi (5 μm)

* For maximum lifetime, recommended maximum pressure for UHPLC particles is 830 bar/12,000 psi.

Intellectual Property: optimizetech.com/patents



Similar to: (5%-phenyl)-methylpolysiloxane

Rxi-5ms Columns (fused silica)

low-polarity phase; Crossbond diphenyl dimethyl polysiloxane

- General-purpose columns for semivolatiles, phenols, amines, residual solvents, drugs of abuse, pesticides, PCB congeners (e.g., Aroclor mixes), solvent impurities.
- Most inert column on the market.
- Tested and guaranteed for ultra-low bleed; improved signal-to-noise ratio for better sensitivity and mass spectral integrity.
- Equivalent to USP G27 and G36 phases.

| ID | df | Length | Temp. Limits | qty. | cat.# |
|----------------|--------------|--------|-------------------|-------|-----------|
| Rxi-5ms | | | | | |
| 0.25 mm | 0.25 μ m | 30 m | -60 to 330/350 °C | ea. | 13423 |
| | 0.25 μ m | 30 m | -60 to 330/350 °C | 6-pk. | 13423-600 |

Topaz GC Inlet Liners

Topaz GC inlet liners feature revolutionary technology and inertness to deliver you the next level of True Blue Performance:

- **Deactivation**—unbelievably low breakdown for accurate and precise low-level GC analyses.
- **Reproducibility**—unbeatable manufacturing controls and QC testing for superior reliability across compound classes.
- **Productivity**—unparalleled cleanliness for maximized GC uptime and lab throughput.
- **100% Satisfaction**—if a liner doesn't perform to your expectations, we will replace it or credit your account.*

Patented

Topaz 4.0 mm ID Single Taper Inlet Liner w/ Wool

for Thermo TRACE 1300/1310 GCs equipped with SSL inlets

| ID x OD x Length | Packing | qty | Similar to Part # | cat.# |
|---|-------------|-------|--------------------------------------|-------|
| Single Taper, Premium Deactivation, Borosilicate Glass | | | | |
| 4.0 mm x 6.5 mm x 78.5 mm | Quartz Wool | 5-pk. | Thermo Fisher Scientific 453A1925-UI | 23447 |

* 100% SATISFACTION GUARANTEE: If your Topaz inlet liner does not perform to your expectations for any reason, simply contact Restek Technical Service or your local Restek representative and provide a sample chromatogram showing the problem. If our GC experts are not able to quickly and completely resolve the issue to your satisfaction, you will be given an account credit or replacement product (same cat.#) along with instructions for returning any unopened product. (Do not return product prior to receiving authorization.) For additional details about Restek's return policy, visit www.restek.com/warranty

California Pesticide Standards

(6 separate mixes)

NEW!

- Meet specific cannabis analysis needs of California set forth by the Bureau of Cannabis Control for regulated category I and II residual pesticide reporting—and of states with similar regulations/programs.
- Ideal for creating multipoint (5-point minimum suggested) calibration curves for GC- and LC-MS/MS.
- Verified composition and stability.
- 66 compounds in 6 x 1 mL ampuls at 100 µg/mL.
- Prepared stock product eliminates the need for in-house standards preparation.
- Restek has the products and expertise you need for cannabis analysis: www.restek.com/cannabis

Each ampul sold separately.

Cat. # 34124: California Pesticide Standard #1 (12 components)

Acephate (30560-19-1)
Chlorpyrifos (2921-88-2)
Coumaphos (56-72-4)
Diazinon (333-41-5)
Dichlorvos (DDVP) (62-73-7)
Dimethoate (60-51-5)
Ethoprophos (13194-48-4)
Malathion (121-75-5)
Methyl parathion (298-00-0)
Mevinphos (7786-34-7)
Naled (300-76-5)
Phosmet (732-11-6)

Cat. # 34125: California Pesticide Standard #2 (11 components)

Abamectin (71751-41-2)
Acequinocyl (57960-19-7)
Bifenthrin (82657-04-3)
Cyfluthrin (68359-37-5)
Cypermethrin (52315-07-8)
Etofenprox (80844-07-1)
Permethrin (*cis* & *trans*) (52645-53-1)
Prallethrin (23031-36-9)
Pyrethrins (8003-34-7)
Spinetoram (J&L) (935545-74-7)
Spinosad (168316-95-8)

Cat. # 34126: California Pesticide Standard #3 (9 components)

Aldicarb (116-06-3)
Bifenazate (149877-41-8)
Carbaryl (Sevin) (63-25-2)
Carbofuran (1563-66-2)
Fenoxycarb (72490-01-8)
Methiocarb (2032-65-7)
Methomyl (16752-77-5)
Oxamyl (23135-22-0)
Propoxur (Baygon) (114-26-1)

Cat. # 34127: California Pesticide Standard #4 (9 components)
Boscalid (188425-85-6)
Captan (133-06-2)
Chlorantraniliprole (500008-45-7)
Daminozide (1596-84-5)
Dimethomorph (110488-70-5)

Fenhexamid (126833-17-8)
Flonicamid (158062-67-0)
Hexythiazox (78587-05-0)
Pyridaben (96489-71-3)

Cat. # 34128: California Pesticide Standard #5 (10 components)

Azoxystrobin (131860-33-8)
Chlorfenapyr (122453-73-0)
Fenpyroximate (111812-58-9)
Kresoxim methyl (143390-89-0)
Metalaxyl (57837-19-1)
Piperonyl butoxide (51-03-6)
Spiromesifen (283594-90-1)
Spirotetramat (203313-25-1)
Spiroxamine (118134-30-8)
Trifloxystrobin (141517-21-7)

Cat. # 34129: California Pesticide Standard #6 (15 components)

Acetamiprid (135410-20-7)
Chlordane (57-74-9)
Clofentazine (74115-24-5)
Etoazoxazole (153233-91-1)
Fipronil (120068-37-3)
Fludioxonil (131341-86-1)
Imazalil (35554-44-0)
Imidacloprid (138261-41-3)
Myclobutanil (88671-89-0)
Paclobutrazol (76738-62-0)
Pentachloronitrobenzene (Quin-tozene) (82-68-8)
Propiconazole (Tilt) (60207-90-1)
Tebuconazole (107534-96-3)
Thiacloprid (111988-49-9)
Thiamethoxam (153719-23-4)

| Conc. in Solvent | CRM? | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | qty. | cat.# |
|---|------|-----------------------------|-----------------------------|---------------------|------------------|------|-------|
| California Pesticide Standard #1 | | | | | | | |
| 100 µg/mL, Acetonitrile, 1mL/ampul | Yes | 6 months | 14 months | Ambient | -20 °C or colder | ea. | 34124 |
| California Pesticide Standard #2 | | | | | | | |
| 100 µg/mL, Acetonitrile, 1mL/ampul | Yes | 6 months | 24 months | Ambient | -20 °C or colder | ea. | 34125 |
| California Pesticide Standard #3 | | | | | | | |
| 100 µg/mL, Acetonitrile, 1mL/ampul | Yes | 6 months | 24 months | Ambient | 0 °C or colder | ea. | 34126 |
| California Pesticide Standard #4 | | | | | | | |
| 100 µg/mL, Acetonitrile, 1mL/ampul | Yes | 6 months | 24 months | Ambient | 0 °C or colder | ea. | 34127 |
| California Pesticide Standard #5 | | | | | | | |
| 100 µg/mL, Acetonitrile, 1mL/ampul | Yes | 6 months | 24 months | Ambient | 0 °C or colder | ea. | 34128 |
| California Pesticide Standard #6 | | | | | | | |
| 100 µg/mL, Acetonitrile, 1mL/ampul | Yes | 6 months | 24 months | Ambient | 0 °C or colder | ea. | 34129 |

Aflatoxins (B1, B2, G1, G2) Standard

(4 components)

NEW!

Ideal for mycotoxin analyses in cannabis and food testing labs.

Aflatoxin B1 (1162-65-8)
Aflatoxin B2 (7220-81-7)

Aflatoxin G1 (1165-39-5)
Aflatoxin G2 (7241-98-7)

| Conc. in Solvent | CRM? | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | Data pack available? | qty. | cat.# |
|------------------------------------|------|-----------------------------|-----------------------------|---------------------|----------------|----------------------|------|-------|
| 10 µg/mL, Acetonitrile, 1 mL/ampul | Yes | 6 months | 36 months | Ambient | 0 °C or colder | No | ea. | 34121 |

Ochratoxin A Standard

Ideal for mycotoxin analyses in cannabis and food testing labs.

NEW!

Ochratoxin A (303-47-9)

| CAS # | Conc. in Solvent | CRM? | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | Data pack available? | qty. | cat.# |
|----------|------------------------------------|------|-----------------------------|-----------------------------|---------------------|----------------|----------------------|------|-------|
| 303-47-9 | 10 µg/mL, Acetonitrile, 1 mL/ampul | Yes | 6 months | 36 months | Ambient | 0 °C or colder | No | ea. | 34122 |

Cannabinoids Standard

(3 components)

U.S. DEA-exempted formulation—no additional customer permits or licensing are required to purchase within the U.S.

Cannabidiol (CBD) (13956-29-1) Cannabinol (CBN) (521-35-7) d9-Tetrahydrocannabinol (d9-THC) (1972-08-3)

| Conc. in Solvent | CRM? | DEA Status | Canadian Test Kit Registration | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | qty. | cat.# |
|---|------|------------|--------------------------------|-----------------------------|-----------------------------|---------------------|-----------------|------|-------|
| 1000 µg/mL each in P&T methanol, 1 mL/ampul | Yes | Exempt | T.K.# 71-048 | 6 months | 24 months | On Ice | 10 °C or colder | ea. | 34014 |

Cannabigerol (CBG) Standard

Excluded from U.S. DEA Controlled Substances Act (CSA) regulatory controls—no customer permits or licensing required to purchase within the U.S.

Cannabigerol (CBG) (25654-31-3)

| Description | CAS # | Conc. in Solvent | CRM? | DEA Status | Canadian Test Kit Registration | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | Type | qty. | cat.# |
|--------------------|------------|--|------|----------------|--------------------------------|-----------------------------|-----------------------------|---------------------|-----------------|--------------|------|-------|
| Cannabigerol (CBG) | 25654-31-3 | 1000 µg/mL in P&T methanol, 1 mL/ampul | Yes | Not Controlled | T.K.# 71-052 | 6 months | 36 months | On Ice | 10 °C or colder | Cannabinoids | ea. | 34067 |

d9-Tetrahydrocannabinol (d9-THC) Standard

U.S. DEA-exempted formulation—no additional customer permits or licensing are required to purchase within the U.S.

d9-tetrahydrocannabinol (d9-THC) (1972-08-3)

| Description | CAS # | Conc. in Solvent | CRM? | DEA Status | Canadian Test Kit Registration | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | Type | qty. | cat.# |
|----------------------------------|-----------|------------------------------------|------|------------|--------------------------------|-----------------------------|-----------------------------|---------------------|-----------------|--------------|------|-------|
| d9-Tetrahydrocannabinol (d9-THC) | 1972-08-3 | 1000 µg/mL in methanol, 1 mL/ampul | Yes | Exempt | T.K.# 71-049 | 6 months | 24 months | On Ice | 10 °C or colder | Cannabinoids | ea. | 34067 |

d9-Tetrahydrocannabinolic Acid A (THCA-A) Standard

U.S. DEA-exempted formulation—no additional customer permits or licensing are required to purchase within the U.S.

d9-tetrahydrocannabinolic acid A (THCA-A) (23978-85-0)

| Description | CAS # | Conc. in Solvent | CRM? | DEA Status | Canadian Test Kit Registration | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | Type | Data pack available? | qty. | cat.# |
|---|------------|--|------|------------|--------------------------------|-----------------------------|-----------------------------|---------------------|------------------|--------------|----------------------|------|-------|
| d9-Tetrahydrocannabinolic acid A (THCA-A) | 23978-85-0 | 1000 µg/mL in acetonitrile, 1 mL/ampul | Yes | Exempt | C.T.K.# 002-011 | 6 months | 21 months | On Ice | -20 °C or colder | Cannabinoids | No | ea. | 34111 |

NEW!

Atrazine-d5 Standard

Isotopically labeled to provide the best approach for pesticide residue quantification.

Atrazine-d5 (163165-75-1)

| Description | CAS # | Conc. in Solvent | CRM? | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | qty. | cat.# |
|-------------|-------------|---------------------------------------|------|-----------------------------|-----------------------------|---------------------|-----------------|------|-------|
| Atrazine-d5 | 163165-75-1 | 100 µg/mL in acetonitrile, 1 mL/ampul | Yes | 6 months | 36 months | Ambient | 10 °C or colder | ea. | 31984 |



Dimethoate-d6 Standard

Isotopically labeled to provide the best approach for pesticide residue quantification.

Dimethoate-d6 (1219794-81-6)

| Description | CAS # | Conc. in Solvent | CRM? | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | qty. | cat.# |
|---------------|--------------|---------------------------------------|------|-----------------------------|-----------------------------|---------------------|-----------------|------|-------|
| Dimethoate-d6 | 1219794-81-6 | 100 µg/mL in acetonitrile, 1 mL/ampul | Yes | 6 months | 36 months | Ambient | 10 °C or colder | ea. | 31988 |

Dichlorvos-d6 Standard

Isotopically labeled to provide the best approach for pesticide residue quantification.

Dichlorvos-d6 (203645-53-8)

| Description | CAS # | Conc. in Solvent | CRM? | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | qty. | cat.# |
|---------------|-------------|----------------------------------|------|-----------------------------|-----------------------------|---------------------|-----------------|------|-------|
| Dichlorvos-d6 | 203645-53-8 | 100 µg/mL in acetone, 1 mL/ampul | Yes | 3 months | 12 months | Ambient | 10 °C or colder | ea. | 31987 |

Carbaryl-d7 Standard

Isotopically labeled to provide the best approach for pesticide residue quantification.

Carbaryl-d7 (362049-56-7)

| Description | CAS # | Conc. in Solvent | CRM? | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | qty. | cat.# |
|-------------|-------------|---------------------------------------|------|-----------------------------|-----------------------------|---------------------|-----------------|------|-------|
| Carbaryl-d7 | 362049-56-7 | 100 µg/mL in acetonitrile, 1 mL/ampul | Yes | 6 months | 31 months | Ambient | 10 °C or colder | ea. | 31985 |

Diazinon-d10 Standard

Isotopically labeled to provide the best approach for pesticide residue quantification.

Diazinon-d10 (diethyl-d10) (100155-47-3)

| Description | CAS # | Conc. in Solvent | CRM? | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | qty. | cat.# |
|----------------------------|-------------|---------------------------------------|------|-----------------------------|-----------------------------|---------------------|-----------------|------|-------|
| Diazinon-d10 (diethyl-d10) | 100155-47-3 | 100 µg/mL in acetonitrile, 1 mL/ampul | Yes | 6 months | 36 months | Ambient | 10 °C or colder | ea. | 31986 |

Diuron-d6 Standard

Isotopically labeled to provide the best approach for pesticide residue quantification.

Diuron-d6 (1007536-67-5)

| Description | CAS # | Conc. in Solvent | CRM? | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | qty. | cat.# |
|-------------|--------------|---------------------------------------|------|-----------------------------|-----------------------------|---------------------|-----------------|------|-------|
| Diuron-d6 | 1007536-67-5 | 100 µg/mL in acetonitrile, 1 mL/ampul | Yes | 6 months | 31 months | Ambient | 10 °C or colder | ea. | 31989 |

Linuron-d6 Standard

Isotopically labeled to provide the best approach for pesticide residue quantification.

Linuron-d6 (1219804-76-8)

| Description | CAS # | Conc. in Solvent | CRM? | Min Shelf Life on Ship Date | Max Shelf Life on Ship Date | Shipping Conditions | Storage Temp. | qty. | cat.# |
|-------------|--------------|---------------------------------------|------|-----------------------------|-----------------------------|---------------------|-----------------|------|-------|
| Linuron-d6 | 1219804-76-8 | 100 µg/mL in acetonitrile, 1 mL/ampul | Yes | 6 months | 31 months | Ambient | 10 °C or colder | ea. | 31990 |

Q-sep QuEChERS Extraction Salts

- Free-flowing salts transfer easily and completely.
- Easy-open packets eliminate the need for a second empty tube for salt transfer.
- Convenient slim packets fit perfectly into tubes to prevent spills.
- Ready-to-use tubes, no glassware required.
- Pre-weighed, ultra-pure extraction salts.
- Ideal for original unbuffered, AOAC (2007.01), and European (EN 15662) QuEChERS methods.

| Description | Material | Method | qty. | cat.# |
|---|---|-------------------|------------|-------|
| Q-sep QuEChERS Extraction Salt Packets Only | 4 g MgSO ₄ , 1 g NaCl, 1 g TSCD, 0.5 g DHS | European EN 15662 | 50 packets | 25849 |

DHS – disodium hydrogen citrate sesquihydrate; MgSO₄ – magnesium sulfate; NaCl – sodium chloride; NaOAc – sodium acetate; TSCD – trisodium citrate dihydrate



25847

Q-sep QuEChERS dSPE Tubes for Extract Cleanup

Fast, Simple Sample Prep for Multiresidue Pesticide Analysis

- Packaged in foil subpacks of 10 for enhanced protection and storage stability.
- Ready-to-use tubes, no glassware required.
- Pre-weighed, ultra-pure sorbents.
- Support original unbuffered, AOAC (2007.01), European (EN 15662), and mini-multiresidue QuEChERS methods.

| Description | Material | Method | Type | Volume | qty. | cat.# |
|---|--|-----------------------------|--|--------|---------|-------|
| Pigmented fruits and vegetables (e.g., strawberries, sweet potatoes, and tomatoes) | | | | | | |
| Q-sep QuEChERS dSPE Tubes | 150 mg MgSO ₄ , 25 mg PSA, 2.5 mg GCB | Mini-multiresidue, EN 15662 | 2 mL Micro-Centrifuge Tubes Prefilled with dSPE Materials for Cleanup (1 mL Extract) | 2 mL | 100-pk. | 26217 |

Note: No entry in the Method column refers to dSPE formulations not specifically included in one of the cited references. These products can be used to accommodate the various needs of specific matrices not directly met by the cited references.



26215

Empty Centrifuge Tubes, Polypropylene

| Description | qty. | cat.# |
|---|--------|-------|
| Empty 50 mL Centrifuge Tube, Polypropylene w/Blue Cap | 50-pk. | 25846 |



25846