

Applications Chart

| ANALYTE CLASS | MATRIX | ANALYTES PER APPLICATION | GRACEPURE™ PRODUCT | PRETREATMENT |
|-------------------------------|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Amphetamines | Urine | Amphetamine and Methamphetamine | C18-Aq, 500mg | Spike urine with 1ug/mL target analytes. Dilute with equal volume of 2% ammonium hydroxide in DI water. |
| Anticonvulsants | Serum | Phenobarbital, Primidone, Carbamazepine, 5,5-Diphenylhydantoin, MPPH (5-Methylphenyl-5-phenylhydantoin) | C18-Low, 500mg | Add 100µL of 0.1M KH ₂ PO ₄ buffer, pH 3.5 to 500µL of serum in a test tube. Add 200µg/mL MPPH, 5-methylphenyl-5-phenylhydantoin as internal standard. Vortex 1 minute. |
| Benzodiazepines | Serum | Norchloriazepoxide, Demoxepam, Chlordiazepoxide, Nitrazepam, Nordiazepam (Metabolite of diazepam), Diazepam | C18-Low, 500mg | Use 500µL serum. Add 500µL internal standard solution: 50µg/mL benzodiazepine. Vortex 1 minute. |
| BHA | Soy Oil | BHA (3-tert-Butyl-4-hydroxyanisole) | Amino, 200mg | Add 10mg BHA into 1mL soy oil and dilute to 10mL with n-Pentane. |
| Caffeine | Coffee | Caffeine | C18-Aq, 500mg | None, will work equally well for any beverage containing caffeine. |
| Carbohydrates | Molasses | Fructose, Glucose, Sucrose | C18-Low, 500mg | Dilute 20g molasses to 250mL with DI water. |
| Carbohydrates | Wine | Ethanol, Glucose, Sucrose | C18-Max, 100mg | None. |
| Chlorinated Pesticides | Water | α-BHC, Lindane, β-BHC, Heptachlor, Aldrin, Heptachlor Epoxide, p,p'-DDE, Dieldrin, o,p'-DDD, Endrin, o,p'-DDT, p,p'-DDD, p,p'-DDT | C18-Fast, 500mg | Due to the large sample volume, attach large volume reservoir to SPE device. |
| Chlorotetracycline | Ointment | Chlorotetracycline | Diol, 500mg | Add 2mL of hexane to 50mg of ointment. Vortex 1 minute. |
| Chlorophenoxy Acid Herbicides | Water | 2,4-D; 2,4,5-T; Silvex | C18-Fast, 500mg | Acidify 100mL water sample to pH 2.2. |
| Desatting | Protein Solution | Cytochrom C, Ribonuclease-A | C18-Aq, 500mg | None. |
| Lactic Acid | Water | Lactic Acid | Anion-X, 500mg | None. |
| Lidocaine, Metabolites | Serum | GX (Glycinexylidide), MEGX (Monoethylglycinexylidide), Lidocaine, Mepivacaine (internal standard) | C18-Low, 500mg | Use 500µL serum. Add 500µL internal standard solution: 50µg/mL Mepivacaine HCl in 0.1M NaH ₂ PO ₄ . Vortex 1 minute. |
| Nitroaromatics and Naphthols | Water | 2,4-DNT, 2-NT, 4-NT, 3-NT, 1-Naphthol, 2-Naphthol | C18-Fast, 500mg | Spike 1000mL tap water with 0.75µg/mL of analytes. |
| Off Flavors | Wine | 4-Ethyl Phenol, 4-Ethyl Gualacol | C18-Low, 500mg | None. |
| Paraben Preservatives | Cosmetics | Methyl Paraben, Propyl Paraben | C18-Low, 500mg | Weigh one gram of cosmetic (hand cream, toothpaste, liquid soap) into a test tube. Add 10mL methanol and vortex one minute. Centrifuge resulting mixture to remove insoluble materials. Remove a 100µL aliquot to a 2mL volumetric flask and dilute to volume with methanol. |
| Perchlorate | Biological Matrix | Perchlorate | Anion-X, 500mg | None. |
| Phenylpropanolamine | Urine | Phenylpropanolamine | C18-Low, 100mg | 1mL urine sample is placed in a small test tube. Add 250mL of carbonate buffer (NaHCO ₃ /Na ₂ CO ₃ , 5:1 w/w) Vortex 1 minute. |
| Phthalate Esters | Drinking Water | Dimethyl Phthalate, Diethyl Phthalate, Diallyl Phthalate, Dibutyl Phthalate, Diamyl Phthalate | C18-Low, 500mg | None. |
| Polyaromatic Hydrocarbons | River Water | Naphthalene, Acenaphthylene, Acenaphthene, Fluorene, Phenanthrene, Anthracene, Fluoranthene, Pyrene, Benz[a]anthracene, Chrysene, Benzo[b]fluoranthene, Benzo[k]fluoranthene, Benzo[a]pyrene, Dibenzo[a,h]anthracene, Benzo[ghi]perylene, Indeno[1,2,3-cd]pyrene | C18-Aq, 500mg | None. |
| Polyaromatic Hydrocarbons | Tap Water | Acenaphthalene, Fluorene, Phenanthrene, Anthracene, Fluoranthene, Pyrene, Benz[a]anthracene, Benzo[b]fluoranthene, Benzo[k]fluoranthene, Benzo[a]pyrene, Dibenzo[a,h]anthracene, Indeno[1,2,3-cd]pyrene, Benzo[ghi]perylene | C18-Low, 500mg | None. |
| Preservatives | Beverages | Propionic Acid, Butyric Acid, Valeric Acid, Caproic Acid, Heptanoic Acid, Caprylic Acid | Anion-X, 500mg | Adjust beverage pH to 10 using KOH. |
| Quinidine | Urine | Quinidine | Diol, 500mg | Add 1mL HCl and 1mL urine sample to a 5mL volumetric flask. Heat to 65°C in a water bath for 10 minutes. Cool and add 1mL ammonium hydroxide. Dilute to volume with distilled water. |
| Salicylic Acid | Urine | Salicylic Acid, Acetylsalicylic Acid | C18-Max, 100mg | Spike 2mL synthetic urine with 100ppm salicylic acid and 100ppm acetylsalicylic acid. |
| Sedatives/Hypnotics | Serum | Barbital, Methypylon, Amobarbital, Phenacetin, Secobarbital, Meprobamate, Glutethimide, Caffeine, Phenobarbital, Methaqualone, Oxazepam, 4-Methyl Primidone, Diazepam, Nodiazepam | C18-Low, 500mg | Use 500µL serum. Add 200µL internal standard solution: 10µg/mL 4-methyl primidone in 0.1M KH ₂ PO ₄ , pH 4. Vortex 1 minute. |
| Steroids | Hydrocortisone Cream | Hydrocortisone | Silica, 500mg | Weigh one gram of cream into a 20mL vial. Add 10mL hexane:ethyl acetate (50:50). Vortex 3 minutes. Decant supernatant into a 50mL volumetric flask. Repeat extraction and combine supernatants. Dilute to volume with hexane:ethyl acetate (50:50). |
| THC | Urine | Δ9-Tetrahydrocannabinol | C18-Low, 500mg | Place 10mL urine sample in a centrifuge tube. Add 0.9mL of 10N NaOH. Cap tube and place in boiling water bath for 15 minutes. Cool to room temperature. Adjust pH to 2. Vortex 1 minute. |
| THC, Metabolites | Urine | Δ9-Tetrahydrocannabinol Methyl Ester, 9-Carboxy-11-nor-Δ9-THC Methyl Ester (Metabolite of #1) | C18-Low, 500mg | Add 1mL methanolic KOH (10% w/v) to 10mL of urine in a test tube. Cap and heat tube to 100°C for 15–20 minutes. Cool to room temperature and adjust pH to 3. |
| Theophylline | Serum | β-Hydroxyethyl Theophylline (internal standard), Theophylline | C18-Low, 100mg | Add 2mL of 0.1M KH ₂ PO ₄ (pH 4) buffer to 1mL serum. Vortex for one minute. |
| Topical Anesthetics | Serum | Benzocaine, Procaine, Mepivacaine | C18-Low, 500mg | Use 500µL serum. Add 500µL internal standard solution: 50µg/mL Mepivacaine HCl in 0.1M NaH ₂ PO ₄ . Vortex 1 minute. |

GracePure™

Applications Chart

| PRECONDITION | LOAD | WASH | ELUTE |
|---------------------------------------------------------|---------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 5mL methanol followed by 5mL DI water. | Apply 10mL sample. | 2mL DI water, followed by 1mL IPA:DI water (25:75). Vacuum 2 minutes. Next wash 1mL hexane, vacuum 2 minutes. Final wash with 1mL IPA. | 3 x 1mL IPA containing 2% ammonium hydroxide. |
| 5mL methanol followed by 5mL DI water. | Add the prepared sample. | 9mL of DI water, vacuum 2 minutes. | 500µL of methanol. |
| 5mL methanol followed by 5mL DI water. | Add the serum sample. | 6mL DI water, vacuum 2 minutes. | 1mL of methanol. |
| 3mL pentane. | Add 1mL sample. | 1.5mL n-pentane. | 2mL ethanol. |
| 5mL methanol followed by 5mL DI water. | Add 1mL prepared sample. | 6mL DI water, vacuum 10 minutes. | 3mL of chloroform. |
| 5mL methanol followed by 5mL DI water. | Add 2mL prepared sample. | No wash, apply vacuum for 5 minutes. | Collect eluate. Filter through a 0.45µm syringe filter. |
| 2mL methanol followed by 2mL DI water. | Add 2mL of wine with the vacuum turned off. | No wash, allow wine to remain in contact with cartridge for 2 minutes. | Turn on vacuum and collect eluant. The organic acids and anthocyanins will remain while the carbohydrates pass through. |
| 5mL methanol followed by 5mL DI water. | Add 100mL of water sample. | No wash, apply vacuum for 5 minutes. | 2mL of ethyl acetate. |
| 3mL of hexane. | Add 500µL prepared sample. | 2mL of hexane, continue vacuum for 3 minutes. | 2mL of a methanol:0.1N HCl solution (50:50). |
| 5mL methanol followed by 5mL DI water. | Add acidified sample. | Wash with 6mL of DI water. | 3mL of chloroform. |
| 3mL methanol followed by 0.025% ammonium hydroxide. | Apply 1mL protein salt solution. | No wash. | 500µL 0.4% TFA followed by 500µL acetonitrile containing 0.4% TFA. Apply vacuum until dry. |
| 2mL 1M NaCl followed by 10mL DI water. | 1mL, 1mL/min. (pH 7). | DI water, 2mL. | 0.1M HCl, 500µL. |
| 5mL methanol followed by 5mL DI water. | Add sample. | 8mL DI water:methanol (75:25), vacuum 2 minutes. | 500µL methanol. |
| 5mL methanol followed by 5mL DI water. | Add 1000mL sample at flow rate of 5mL/min. | No wash. | Elute with 3 x 1mL methanol:water (50:50). Air dry after each elution. |
| 5mL methanol followed by 5mL DI water. | Apply 10mL wine sample. | 5mL water. | 1mL isopropyl alcohol. |
| 5mL methanol followed by 5mL DI water. | Add 2mL prepared sample. | 3mL DI water, vacuum 2 minutes. | 1mL methanol. |
| 3mL 0.5M NaCl followed by 3mL DI water. | Apply 1mL sample. | No wash. | 3 x 0.75mL of 0.1M NaCl. |
| 2mL methanol followed by 2mL DI water. | Add the buffered urine. | 2mL DI water, vacuum 2 minutes. | 6mL of chloroform:isopropanol (90:10) through the cartridge. Repeat with an additional 0.2mL. |
| 5mL methanol followed by 5mL DI water. | Add 200mL water sample. | 3mL DI water. | Pass two 500µL aliquots of ethyl acetate. |
| 5mL methanol followed by 5mL DI water. | Apply 200mL water containing PAH's. | 2mL DI water followed by 2mL IPA:Water (20:80). | 2 x 2mL methanol. |
| 6mL 2-propanol:DI water (15:85). | Add 100mL water sample. | 2mL 2-propanol:DI water (15:85). | 1mL methylene chloride. |
| 10mL DI water. | Apply 8mL beverage sample. | 20mL DI water. | 1mL 1.0N HCl followed by 1mL methanol. |
| 3mL methanol followed by 3mL DI water adjusted to pH 9. | Add 500µL prepared sample. | 1mL of distilled water, continue vacuum for 2 minutes to remove residual wash solution. | Pass two aliquots of 500µL methanol. |
| 3mL methanol followed by 3mL DI water. | Add 2mL spike urine. | 2mL 50mM phosphate buffer monobasic, pH 2. | 2mL methanol:water (50:50). |
| 5mL methanol followed by 5mL DI water. | Add prepared serum sample. | 6mL DI water, vacuum 2 minutes. | 500µL acetone. |
| 2mL, hexane:acetone (80:20). | Add 1mL prepared sample. | 2mL of hexane:acetone (80:20) vacuum 2 minutes. | Pass two aliquots of 500µL methanol. |
| 5mL methanol followed by 5mL DI water. | Add prepared urine sample. | Wash first: 10mL of 0.1M HCl. Wash second 25mL of 50µM phosphoric acid containing 10% acetonitrile. Vacuum 2 minutes. | 3mL of acetone through the cartridge. Collect eluate and add 1.5mL of methylene chloride, centrifuge 5 minutes. Remove upper phase and add 1.5mL of hexane. Centrifuge for 5 minutes. Remove upper phase once again and dry the treated sample. Redissolve in 200µL of chloroform for subsequent GC analysis. |
| 5mL methanol followed by 5mL DI water. | Add prepared urine sample. | 5mL DI water followed by 5mL of acetonitrile:water (40:60). Vacuum 2 minutes. | 2mL of methanol. |
| 2mL methanol followed by 2mL DI water. | Add buffered serum. | 2mL DI water, vacuum 2 minutes. | 1mL of methanol. |
| 5mL methanol followed by 5mL DI water. | Add sample. | 8mL DI water:methanol (75:25), vacuum 2 minutes. | Pass 500µL of methanol and dry. Redissolve in 200µL of chloroform for subsequent analysis by gas chromatography. |

GracePure™