

**Abstract**


C<sub>18</sub> is the most popular stationary phase used in reversed-phase chromatography, but has difficulty in retain highly polar compounds. Capcell Pak ADME (ADME) is developed by firstly introducing adamantyl groups as the functional group of the reversed stationary phase of a HPLC column. The changed alkyl structure of adamantyle group provides ADME separation characteristics with a balance between moderate hydrophobicity and extra-high surface polarity which is completely different from the conventional C<sub>18</sub>. In this study, the characteristics of ADME and its applications compared with C<sub>18</sub> columns in the analysis of metabolite are reported.

Key words: Adamantyl functional Group, polar compound, metabolite, HPLC, LC-MS.

**Experiment & Results**

**Characteristics of ADME**

The hydrophobicity and Surface Polarity of Capcell Pak ADME were measured with 10 standard-reagents method (S. Kobayashi etc. *J. Chromatogr. A*, **1998**, 828, 75-81) and plotted in the parameter graph among Capcell Pak series, Capcell Core series and other conventional C<sub>18</sub> columns (Fig.1). Capcell Pak ADME positions at extra-high surface polarity and moderate hydrophobicity which is unavailable with any other polar columns and conventional C<sub>18</sub> columns. Based on the Capcell Pak technology, the implementation rate is tightly controlled at an optimized density in order to make the performance of ADME group in maximum.



[Fig. 1a] Parameter of hydrophobicity and polarity


**Physical characteristics of ADME (Table 1)**

Particle Size (μm)	Micro pore Diameter (μmm)	Special Surface (m <sup>2</sup> /g)	Density (μmol/m <sup>2</sup> )	C %	pH range
5	10	300	2.8	12	2-9
3	10	300	2.8	12	2-9


**Capcell Pak ADME Vs C18 columns**

With the novel balance between hydrophobicity and polarity never experienced by other reverse-phase columns, ADME provide improved separation to high polar compounds such as the metabolites. In the simultaneous analysis of metabolites of caffeine (Fig.2) and tryptophan (Fig.3), ADME provided improved separation compared with the conventional high polar columns and C<sub>18</sub> columns.

**Improved separation with ADME (isocratic)**



**Improvement in retention, Rs, and NTP**





[Fig. 2] Chromatograms of Xanthenes

## Capcell Pak ADME Vs C<sub>18</sub> columns

Capcell Pak ADME strongly retains the highly polar compounds in the process of metabolism of the tryptophan (Fig 3a) due to its large surface polarity, and provides improved separation even to the compounds hard to retain by the conventional C<sub>18</sub> columns. Especially, ADME could advance the retention of high polar compounds much more under a gradient condition. ADME provided excellent separation capacity to retain the highly polar metabolites of tryptophan with a wide analysis window. In addition, downsizing of particle size to 3μm further improved the number of theory plate (NTP) and resolution.(Fig. 3b).



Improved separation of metabolites ( gradient)



[Fig. 3a] Metabolism pathway of tryptophan

[HPLC Conditions]	
Column size	: [S3] 4.6 mm I.d. x 100 mm, [S5] 4.6 mm I.d. x 150 mm
Mobile phase	: A) 10 mmol/L Ammonium formate buffer (pH4) B) CH <sub>3</sub> CN
[S3] B	: 5% (0 min) → 25% (13.3 min) → 5% (13.4 min)
Gradient	
[S5] B	: 5% (0 min) → 25% (20.0 min) → 5% (20.1 min)
Gradient	
Flow rate	: 1 mL/min
Temperature	: 40 °C
Detection	: UV 254 nm
Inj. vol.	: 1 μL
Sample	: 1. Uracil, 2. 3-Hydroxykynurenine, 3. Kynurenine, 4. Tryptophan, 5. 3-Hydroxyanthranilic acid, 6. Anthranilic acid


[Fig. 3b] Chromatograms of tryptophan and its metabolites

## Lineup of Capcell Core ADME

9 steroid hormone in the Cholesterol biosynthetic pathway were certainly separated with Capcell Pak ADME (5μm) and its core-shell type. Fast simultaneous analysis was reached simply by switching ADME to Capcell Core ADME (2.7μm) with arrows 2-fold linear velocity and short gradient program. (Fig. 4)

[HPLC Conditions]	
Column	: ① CAPCELL PAK ADME S5; 4.6 mm I.d. x 150 mm ② CAPCELL CORE ADME S2.7; 2.1 mm I.d. x 50 mm
Mobile phase	: A) 5 mmol/L HCOONH <sub>4</sub> B) CH <sub>3</sub> CN
	: ① B 20% (0 min) → 55% (30 min) → 20% (30.1 min) Gradient ② B 20% (0 min) → 55% (10 min) → 20% (10.1 min) Gradient
Flow rate	: ① 1 mL/min, ② 400 μL/min
Temperature	: 40 °C; Detection : UV 220 nm
Inj. vol.	: ① 5 μL, ② 25 μL
Sample	: 1. Aldosterone, 2. Cortisol, 3. Cortisone, 4. Corticosterone, 5. 17β-Estradiol, 6. Testosterone, 7. Estrone, 8. Ethinylestradiol, 9. 17α-Hydroxyprogesterone

[Fig. 4] Chromatograms of tryptophan and its metabolites



## Conclusion & Discussion

Capcell Pak ADME and Capcell Core ADME firstly introduced the adamantyl group as the functional group of reverse-phase columns and provide improved separation of highly-polar compounds. ADME is expected as an alternative to C<sub>18</sub> column in simultaneous analysis of samples having diverse polarities such as metabolites and fast separation of polar compounds.

Website & contact

Frontier Science Business division  
Shiseido Co., Ltd.

1-6-2 Higashi-shimbashi, Tokyo, Japan  
Email: [hplc@to.shiseido.co.jp](mailto:hplc@to.shiseido.co.jp)

URL: <http://hplc.shiseido.co.jp/e>

SHISEIDO