



OSAKA SODA

CAPCELL PAK

ADME-HR / INERT ADME-HR



 OSAKA SODA CO., LTD.

A Cage-Structured C₁₂ column Achieving Strong Retention for Polar Analytes under Reversed-Phase Mode

Best second choice

CAPCELL PAK ADME-HR CAPCELL PAK INERT ADME-HR

The introduction of Adamantylethyl groups provides a hydrophobic interaction while maintaining high surface polarity, resulting in retaining polar analytes even under water-rich mobile phases.

Physical property values

CAPCELL PAK ADME-HR

Particle Size (μm)	Pore Size (nm)	Surface Area (m ² /g)	Ligand Density (μmol/m ²)	C%	pH Range	Max. Pressure (MPa)
2	10	310	2.7	12	2~9	100
3	10	310	2.7	12	2~9	20
5	10	310	2.7	12	2~9	20

CAPCELL PAK INERT ADME-HR

Particle Size (μm)	Pore Size (nm)	Surface Area (m ² /g)	Ligand Density (μmol/m ²)	C%	pH Range	Max. Pressure (MPa)
3	10	310	2.7	12	2~9	50

Structure of INERT Column




What is ADME?

ADME is an abbreviation of **Adamantane**, which consists of ten carbons in a diamond-like structure. Ethyl groups are introduced to the **Adamantane** as a spacer and employed as a unique bonded phase for the CAPCELL PAK ADME-HR columns.



Adamantylethyl Groups (ADME Groups)

Parameter Map





A Truly Unprecedented Balance of Hydrophobicity and Surface Polarity


Comparison of Stereoselectivity

HPLC Conditions

- Column size : 2.1 mm i.d. x 50 mm
- Mobile phase: H₂O / CH₃CN = 70 / 30
- Flow rate : 0.4 mL/min
- Temperature : 40 °C
- Detection : PDA 220 nm
- Inj. vol. : 3 μL (50 μg/mL each)
- Sample : 1. 17β-Estradiol (E2)
2. Estrone (E1)
3. Ethynylestradiol (EE)



Comparison of Surface Polarity between a C₁₈ and ADME column



The caged-structured ADME groups offer interaction with the surface of the silica, resulting in providing a unique selectivity compared to a conventional C₁₈ column.


Highly Stable under 100 % Water Mobile Phases

As shown on the right, the efficiency is very stable even after 600 injections under acidic mobile phase.


HPLC Conditions

Column size : 2.1 mm i.d. x 150 mm
 Mobile phase: 10 mmol/L HCOONH₄ (Adjusted with formic acid at pH 3)
 Flow rate : 0.2 mL/min
 Temperature: 40 °C
 Detection : UV 254 nm
 Inj. vol. : 1 µL
 Sample : 1. Cytosine 2. Uracil 3. Guanine 4. Hypoxanthine
 5. Xanthine 6. Oxipurinol 7. Allopurinol 8. Thymine

Results after 600 Consecutive Injections



Variation of Efficiency of Thymine




Retention Behavior of Biogenic Amines under 100 % Water Mobile Phase

As shown below, CAPCELL PAK ADME-HR offers stronger retention of polar analytes, resulting in delivering complete separation for all analytes.


CAPCELL PAK ADME-HR S3
13.6 MPa




CAPCELL PAK ADME-HR S5
8.1 MPa



Highly Polar C₁₈ 5 µm column
Pressure: 7.2 MPa



Hybrid ODS 5 µm column
Pressure: 7.9 MPa



$N_6=20967$


$N_6=16031$

$N_4=12377$

$N_6=12460$

HPLC Conditions



Column size : 4.6 mm i.d. x 150 mm
 Mobile phase: 0.1 vol% HCOOH
 Flow rate : 1.0 mL/min
 Temperature : 40 °C
 Detector : NQAD (Evaporation 60 °C, Nebulizer 30 °C)
 Inj. vol. : 3 µL





Analysis of Nucleic-Acid Bases and Nucleosides via Gradient Elution

As shown below, the unique selectivity provided from CAPCELL PAK ADME-HR shows complete separation for all analytes, but with stronger retention of polar analytes.


CAPCELL PAK ADME-HR S5



Highly Polar C₁₈ 5 µm column



Hybrid ODS 5 µm column



HPLC Conditions

Column size : 2.0 or 2.1 mm i.d. x 150 mm
 Mobile phase: A) 10 mmol/L HCOONH₄, H₂O
 B) CH₃CN
 B 1 % (0 min) -> 1 % (1 min) ->
 40 % (15 min) -> 1 % (15.1 min) Gradient
 Flow rate : 0.2 mL/min
 Temperature : 40 °C
 Detection : UV 254 nm
 Inj. vol. : 2 µL

Ordering Information

CAPCELL PAK ADME-HR

P/N	Description	Particle Size (μm)	I.D. (mm)	Length (mm)	List Price (JPY)
93300	ADME-HR	2	2.1	20	*
93301	ADME-HR	2	2.1	50	*
93302	ADME-HR	2	2.1	100	*
93310	ADME-HR(1/32)	3	0.3	100	110,000
93311	ADME-HR(1/32)	3	0.3	150	120,000
93312	ADME-HR	3	1.0	100	76,000
93320	ADME-HR	3	2.1	20	59,000
93321	ADME-HR	3	2.1	35	54,000
93322	ADME-HR	3	2.1	50	56,000
93323	ADME-HR	3	2.1	75	61,000
93324	ADME-HR	3	2.1	100	65,000
93325	ADME-HR	3	2.1	150	72,000
93326	ADME-HR	3	2.1	250	88,000
93330	ADME-HR	3	3.0	50	54,000
93331	ADME-HR	3	3.0	100	63,000
93332	ADME-HR	3	3.0	150	77,000
93340	ADME-HR	3	4.6	35	47,000
93341	ADME-HR	3	4.6	50	52,000
93342	ADME-HR	3	4.6	75	56,000
93343	ADME-HR	3	4.6	100	60,000
93344	ADME-HR	3	4.6	150	69,000
93345	ADME-HR	3	4.6	250	83,000
12600	ADME-HR CARTRIDGE (2PCS)	3	2.0	10	30,000
12601	ADME-HR CARTRIDGE (2PCS)	3	4.0	10	30,000
12415	CARTRIDGE HOLDER 10 (L)	-	-	10	20,000

P/N	Description	Particle Size (μm)	I.D. (mm)	Length (mm)	List Price (JPY)
93350	ADME-HR	5	2.1	20	59,000
93351	ADME-HR	5	2.1	35	54,000
93352	ADME-HR	5	2.1	50	56,000
93353	ADME-HR	5	2.1	75	61,000
93354	ADME-HR	5	2.1	100	65,000
93355	ADME-HR	5	2.1	150	72,000
93356	ADME-HR	5	2.1	250	88,000
93360	ADME-HR	5	3.0	150	77,000
93361	ADME-HR	5	3.0	250	99,000
93370	ADME-HR	5	4.6	35	47,000
93371	ADME-HR	5	4.6	50	52,000
93372	ADME-HR	5	4.6	75	56,000
93373	ADME-HR	5	4.6	100	60,000
93374	ADME-HR	5	4.6	150	69,000
93375	ADME-HR	5	4.6	250	83,000
93380	ADME-HR	5	10	35	100,000
93381	ADME-HR	5	10	150	210,000
93382	ADME-HR	5	10	250	270,000
93390	ADME-HR	5	20	35	170,000
93391	ADME-HR	5	20	50	240,000
93392	ADME-HR	5	20	100	340,000
93393	ADME-HR	5	20	150	450,000
93394	ADME-HR	5	20	250	750,000
12610	ADME-HR CARTRIDGE (2PCS)	5	2.0	10	30,000
12611	ADME-HR CARTRIDGE (2PCS)	5	4.0	10	30,000
12415	CARTRIDGE HOLDER 10 (L)	-	-	10	20,000

CAPCELL PAK INERT ADME-HR

P/N	Description	Particle Size (μm)	I.D. (mm)	Length (mm)	List Price (JPY)
95001	ADME-HR	3	2.0	50	74,000
95002	ADME-HR	3	2.0	100	83,000
95003	ADME-HR	3	2.0	150	90,000

* 2 μm columns will be available by the end of 2018

CONTACT US :

OSAKA SODA CO., LTD.

DAISO Fine Chem USA, Inc.

DAISO Fine Chem GmbH

1-12-18 Awaza, Nishi-ku, Osaka 550-0011, Japan
TEL: +81-(0)6-6110-1598, FAX: +81-(0)6-6110-1612
E-mail: silica@osaka-soda.co.jp

3858 Carson Street, Suite 126, Torrance, CA 90503, USA
TEL: +1-310-540-5312, FAX: +1-310-540-5332
E-mail: daiso_us@osaka-soda.co.jp

Immermannstrasse 13, 40210, Düsseldorf, Germany
TEL: +49-(0)211-83025168, FAX: +49-(0)211-83025213
E-mail: gmbh@osaka-soda.co.jp

© Copyright OSAKA SODA CO., LTD. March 2006. All rights reserved.
No part of this catalog may be copied or reprinted, for commercial purposes, without prior written permission from OSAKA SODA CO., LTD.