

For Protein Separation

Proteonavi

SHISEIDO

Contents

1	Protein separation in reversed-phase mode	...1
2	Properties of Proteonavi	...2
2-1	Property values	...2
2-2	Evaluation of Hydrophobicity and surface polarity	...2
3	Features of Proteonavi	...4
3-1	Specifically large retention of proteins	...4
3-2	Outstanding acidic durability	...6
3-3	Minimal protein adsorption and high column efficiency	...8
3-4	Excellent lot reproducibility	...9
3-5	Minimal column bleed	...10
4	Scale up from analytical to preparative size	...11
4-1	Scale up process	...11
4-2	Separation efficiency, column size, and sample load	...14
5	Application Data	...16
5-1	Insulin	...16
5-2	Thiogloblin	...17
6	Product List	...18
7	Global distributors	...19

1 Protein separation in reversed-phase mode

Chromatographic separations of proteins are based on their difference in either, size (size exclusion chromatography), electric charges (ion exchange chromatography), or hydrophobicity (reversed-phase chromatography).

“Proteonavi” was developed to separate proteins and peptides in reversed phase mode. In reversed-phase mode, it is known that proteins and peptides with higher-order structure are denatured to some extent in course of being retained on the stationary phase after introduction to the column. Once the retention with denaturation is occurred, molecules will not migrate until the organic content of the mobile phase is raised to an appropriate level under a gradient program.

In reversed-phase separation, peak shapes are often deteriorated when secondary interactions, such as electrostatic interactions with partial bare silica or metal impurities, exist between the stationary phase and proteins. Listed below are common facts in analytical and preparative separations of proteins and peptides.

1. Stationary phase with long alkyl chains (C18) shows irreversible adsorption of protein.
2. Stationary phase with shorter alkyl chains (C1, C4, or C8) have a lower durability under acidic mobile phases.
3. Stationary phase with shorter alkyl chains are influenced by silica, as a starting material, showing a lot variation.
4. Stationary phase with shorter alkyl chains often shows insufficient separation.

All of the above issues are related to chemistry to introduce alkyl chains to silica. The type of alkyl chain and the way to anchor it on silica seem to be the keys of preparing an appropriate stationary phase for protein and peptide analysis. Proteonavi has a C4-modified surface designed after comprehensive basic research.



Shiseido takes pride in introducing Proteonavi for all researchers!

2 Properties of Proteonavi

2-1 Property values

Proteonavi utilizes high-purity silica with few metal impurities, and shows minimal irreversible adsorption for proteins and peptides. Its pore size is as wide as 30 nm, so that large proteins are able to have enough interactions with the stationary phase.


Property values

Functional group	Pore size (nm)	Particle size (μm)	Specific surface area (m^2/g)	Carbon content (%)	Alkyl group density ($\mu\text{mol}/\text{m}^2$)	pH	USP
C ₄	30	5	105	3.3	3.2	2 - 10	L26

2-2 Evaluation of hydrophobicity and surface polarity

Nature of reversed-phase columns can be studied by measuring separation factors of different pairs of standard compounds. The strength of hydrophobicity can be expressed as a separation factor of toluene and benzene, showing a retention increase generated by one methyl group. Similarly, surface polarity is calculated between methyl benzoate and benzene, observing an influence of ester moiety on retention.

Results obtained with CAPCELL PAK SG 300 C₁, C₈, and C₁₈, Proteonavi and two other commercially available columns for protein separation were plotted on the hydrophobicity/surface polarity graph below.



Hydrophobicity is the most fundamental parameter in reversed-phase columns. Generally, a stationary phase with a large carbon content (larger alkyl density and/or longer alkyl chain) shows a large value. A large hydrophobicity is preferred when a

large retention or a high loadability is required. A large value, however, may not be always good for macromolecules with a structure of higher degrees, such as proteins, since a long alkyl chain (C18) is more likely to cause an irreversible adsorption accompanied by their denaturation.

Surface polarity contributes to the retention and separation of polar compounds, and therefore, generates selectivity. Residual silanol groups on silica surface are thought to be one of the sources of surface polarity. While peak shapes of neutral compounds are hardly influenced by them, acidic silanol groups, if any, are known to cause a peak tailing for those having a basic moiety, such as proteins.

CAPCELL PAK SG300 columns form a characteristic curve in the graph, showing the influences of alkyl chain length on both values. However, even C1, the one having the shortest chain, shows smaller polarity values than those of other products (V and Y), which can be attributed to the polymer-coated surface of CAPCELL PAK series.


Protonavi shows the surface polarity close to that of CAPCELL PAK C1, in spite of its C4 structure. Its entirely different surface coating process provides a new balance of hydrophobicity and surface polarity for protein and peptide analysis.

3. Features of Proteonavi

3-1 Specifically large retention of proteins

Proteonavi, having C4 chains, showed overall retention of standard neutral small compounds, which is somewhat between C1 and C8 columns, as expected. However, it retains six standard proteins more than the C8 column, apparently because of its structure designed specifically for protein separation. The large retention does not contradict with the unique property values discussed at the section before.


Neutral low molecular organic compounds




【HPLC Conditions】

Mobile phase : CH₃CN / H₂O = 40 / 60
Column size : 2.0 mm i.d. x 150 mm
Flow rate : 200 μL/min
Temp. : 40 °C
Detection : UV 254 nm
Inj. vol. : 2 μL

Standard proteins



Standard peptides



Standard proteins

1. Ribonuclease A (M.W. : 13.7 kDa)
2. Cytochrome C (M.W. : 12.4 kDa)
3. Lysozyme (M.W. : 14.3 kDa)
4. B.S. albumin (M.W. : 69.0 kDa)
5. Myoglobin (M.W. : 17.0 kDa)
6. Ovalbumin (M.W. : 45.0 kDa)

【HPLC Conditions】

Column size : 2.0 mm i.d. x 150 mm
 Mobile phase : A) 0.1vol% TFA / H₂O
 B) 0.1vol% TFA / CH₃CN
 B 20% (0 min) → 70% (25 min) Gra.
 Flow rate : 200 μL/min
 Temp. : 40 °C
 Detection : UV 220 nm
 Inj. vol. : 2 μL
 Sample : 6 Proteins (STD)
 Equal amount of 1000μg/mL for 1 to 3 and
 5, 2000μg/mL for 4, and 3000μg/mL for 6
 were mixed.

Standard peptides

1. Methionine-Enkephalin (M.W. 574)
2. Bradykinin (M.W. 1060)
3. Leucine-Enkephalin (M.W. 556)
4. Angiotensin II (M.W. 1046)
5. Neurotensin (M.W. 1673)
6. Angiotensin I (M.W. 1296)

【HPLC Conditions】


Column size : 2.0 mm i.d. x 150 mm
 Mobile phase : A) 0.1vol% TFA / H₂O
 B) 0.1vol% TFA / CH₃CN
 B 15% (0 min) → 30% (20 min) Gra.
 Flow rate : 200 μL/min
 Temp. : 40 °C
 Detection : UV 220 nm
 Inj. vol. : 2 μL
 Sample : 6 Peptides (STD) 100μg/mL each

3-2 Outstanding acidic durability

Acidic mobile phases containing trifluoroacetic acid (TFA) are often used in protein separation in reversed-phase chromatography. Acidic hydrolysis of alkyl chain is known to occur more or less on any silica-based columns, which is more pronouncing for those with shorter alkyl chains. Partial loss of alkyl groups generally causes a decrease of overall retention of analytes and a change of selectivity toward them.

In order to evaluate durability of different columns under acidic conditions, an accelerated test was performed in the conditions shown below. Retention time change of benzylalcohol was examined under the strongly acidic mobile phase at the elevated temperature.


As shown in the graph, all of the columns tested showed loss of retention as time passed. While the column Y showed a rapid drop, and V lost approximately 30% of its retention during the process, Proteonavi kept more than 90% of the initial retention time until the end.



【HPLC Conditions】

Column size : 4.6 mm i.d. x 150 mm
 Mobile phase : 1 vol%TFA, CH₃OH / H₂O = 10 / 90
 Flow rate : 1 mL/min
 Temp. : 60 °C
 Detection : UV 254 nm
 Inj. vol. : 7μL
 Sample : Uracil (50 μg/mL), Benzylalcohol (1000 μg/mL)
 Runtime : 60 min per run

Actual chromatograms of the first and 20th runs (10th run for the column Y) are shown below. The absence of change in uracil peak (shape and position) at t_0 indicates that the test hardly damaged silica support in all of the columns. Retention time of benzyl alcohol, on the contrary, showed a significant difference from one column to another.




3-3 Minimal protein adsorption and high column efficiency

The primary source of protein retention on reversed phase is hydrophobic interaction between hydrophobic parts of amino acids residues and alkyl chains in the stationary phase. When secondary interactions, such as coulombic interactions by acidic silanols or metal impurities in the stationary phase, exist to a large extent, peak shapes of proteins will be deteriorated due to their slow kinetics. Proteonavi shows better peak profiles and higher resolution among standard proteins, in comparison with conventional wide-pore columns.

【HPLC Conditions】

Column size : 2.0 mm i.d. x 150 mm
 Mobile phase : A) 0.1 vol% TFA / H₂O
 B) 0.1 vol% TFA / CH₃CN
 B 20% (0 min) → 70% (15 min) Gra.
 Flow rate : 200 μL/min
 Temp. : 40 °C
 Detection : UV 220 nm
 Inj. vol. : 2 μL
 Sample : 6 Proteins (STD)

1. Ribonuclease A (M.W. : 13.7 kDa)
 2. Cytochrome C (M.W. : 12.4 kDa)
 3. Lysozyme (M.W. : 14.3 kDa)
 4. B.S. albumin (M.W. : 69.0 kDa)
 5. Myoglobin (M.W. : 17.0 kDa)
 6. Ovalbumin (M.W. : 45.0 kDa)



3-4 Excellent lot reproducibility

Lot-to-lot variation in separation is often discussed in protein analysis. Proteonavi has excellent reproducibility of retention time and separation.

【HPLC Conditions】

Column size : 2.0 mm i.d. x 150 mm
Mobile phase : A) 0.1 vol% TFA / H₂O
 B) 0.1 vol% TFA / CH₃CN
 B 20% (0 min) → 70% (15 min) Gradient

Flow rate : 200 μL /min


Temp. : 40 °C

Detection : UV 220 nm

Inj. vol. : 2 μL

Sample : 6 Proteins (STD)


- | | |
|-------------------|-------------------|
| 1. Ribonuclease A | (M.W. : 13.7 kDa) |
| 2. Cytochrome C | (M.W. : 12.4 kDa) |
| 3. Lysozyme | (M.W. : 14.3 kDa) |
| 4. B.S. albumin | (M.W. : 69.0 kDa) |
| 5. Myoglobin | (M.W. : 17.0 kDa) |
| 6. Ovalbumin | (M.W. : 45.0 kDa) |




3-5 Minimal column bleed

“Column bleed” used to be a term used in gas chromatography with universal detectors, such as flame ionization detector. The recent development of interfacing technologies in mass spectrometry brought the term into the field of LC. Column bleed of an LC column is a synonym of a signal level of total ion chromatogram (TIC) obtained with LC-MS under gradient conditions at full scan mode.

The source of column bleed is a bonded phase on silica. It may affect ionization efficiency of analytes of interest, or will be potential impurities in preparative separation. Proteonavi was designed to minimize the level of column bleed. A comparison among conventional wide-pore columns and Proteonavi is shown below. These background chromatograms were obtained under gradient conditions without sample introduction.

**【HPLC conditions】**

Column size : 2.0 mm i.d. x 150 mm
Mobile phase : A) 0.1 vol% TFA / H₂O,
B) 0.1 vol% TFA / CH₃CN
Flow rate : 200 μL/min
Detection : LC-MS (AccuTOF)



4 Scale up from analytical to preparative size

4-1 Scale up process

Proteonavi is suitable to large-scale preparative separation, taking advantage of its low-column bleed nature and high separation efficiency.

The followings are a common process to set up conditions for preparation separation.

Step
1

Pre-optimizing conditions in analytical scale

- 1) Mobile phase optimization
- 2) Know temperature dependence of separation
- 3) Study loadability by increasing sample amount.

Step
2

Determination of size

Select the inner diameter and particle size of a preparative column.

Step
3

Optimizing conditions in preparative scale

- 1) Study sample load and peak resolution
- 2) Column pressure and flow rate

Step
4

Running real samples

The table below gives a conversion factor to calculate flow rate or sample size for each column size under the same linear velocity. In practice, a flow rate tends to be lowered in a large scale separation because of a pump capacity or safety reasons.


Conversion factors for flow rate and sample size

Inner diameter of column (mm i.d.)	Relative Flow rate or sample size
4.6	1
10	5
20	19
30	43

Here are the tips in preparative separation.


i) Impurities on the edges

When impurities exist on the either edge of target peak, sample load and timing of fractionation will be determined by compromise between throughput and purity.




ii) Purity

Purity of collected sample is always inversely related to recovery or throughput.





iii) Throughput

**How to obtain the maximum production efficiency?**

Although quite a few theoretical works have been published in preparative separation, optimization of fractionation with purity evaluation of each fraction is always necessary. A composition of mixture and target purity differs in every case. Theory does not help you much (-_-).

4-2 Separation efficiency, column size, and sample load


Separation efficiency at the same sample load was compared using Proteonavi columns of 150-mm and 250-mm long. The 250-mm column shows a higher resolution when an appropriate gradient time program is used. A resolution between peaks tends to be lowered as a sample size increases. However, the loss of separation efficiency by increased sample load will be moderate in a longer column.



【HPLC Conditions】

Mobile phase : A) 0.1 vol% TFA / H₂O, B) 0.1 vol% TFA / CH₃CN
 Flow rate : 1 mL/min
 Temperature : 40 °C
 Detection : 220 nm (Diodearray)
 Sample dissolved in : H₂O

Separation of two proteins of various sample loads were compared. (Detection wavelength was set at 225nm in order to avoid optical saturation in detection) The two proteins were baseline-separated even at 20-fold sample load in this particular case.



【HPLC Conditions】


Column : Proteonavi , 4.6 mm i.d. x 150 mm
 Mobile phase : A) 0.1 vol% TFA / H₂O, B) 1vol% TFA / CH₃CN
 B% 20% (0 min) → 70% (15 min) → 70% (16 min)
 → 20% (16.1 min) Gra.
 Flow rate : 200 µL/min
 Temp. : 40 °C
 Detection : 225 nm (Diodearray)
 Inj. vol : 10 µL, 50µL, 100µL, 200µL (1000 µg/mL each)
 Sample : Bovine serum albumin (BSA), Myoglobin
 Sample dissolved in : H₂O

5 Application Data (FAQ)

5-1 Insulin (5.7 kDa)



Q. Can we analyze insulin with Proteonavi?

A. Shown below are the analytical results using Proteonavi. Insulin was detected by photodiodearray (PDA) and mass spectroscopy (MS).




【HPLC Conditions】

Column : Proteonavi ; 2.0 mm i.d. x 150 mm
 Mobile phase : A) 0.1 vol% TFA / H₂O
 B) 0.1 vol% TFA / CH₃CN
 B% 20% (0min) → 50% (15min) →
 20% (15.1min) → 20% (25min) Gra.
 Flow rate : 200 μL/min
 Temp. : 40 °C
 Detection : 224 nm (Diodearray)
 Inj. vol. : 2 μL (100 μg/mL)
 Sample dissolved in : M. phase



Protein, which is constituted by amino acids, is detected as multi-charged ions.



Isotope pattern

Particular pattern of large molecule by ¹³C is obtained.
 (Resolution: approximately 100,000)


【LC-MS Conditions】

Column : Proteonavi ; 2.0 mm i.d. x 150 mm
 M.phase : A) 0.1 vol% HCOOH / CH₃CN = 95 / 5
 B) 0.1 vol% HCOOH / CH₃CN = 20 / 80
 B% 10% (0min) → 70% (10min) →
 10% (10.1min) → 10% (20min) Gra.
 Flow rate : 200 μL/min
 Temperature : 40 °C
 Detection : ESI, Positive (LTQ-Orbitrap)
 Inj. amount : 1 nmol
 Sample dissolved in : Eluent

5-2 Thiogloblin (660kDa)

Q. How large molecules can Proteonavi separate?

A. There is no clear maximum molecular weight to be generally applied to all the proteins. However, very large molecules, such as thiogloblin, may not show a sharp peak profile any more. In addition, heterogeneity of the protein leads to the characteristic multi-peak profile in a chromatogram. Nevertheless, Proteonavi shows a slightly larger retention time and better separation between the last two major peak groups among the tested.



6 Product List

■ Analytical size

Product number	Particle size (μm)	Inner diameter (mm)	Length (mm)
80200	5	2.0	50
80201			150
80202			250
80203		4.6	50
80204			150
80205			250
11530		2.0	10 (cartridge, 2pcs)
11531		4.0	10 (cartridge, 2pcs)
12415		—	—

※ Cartridge holder is needed when using cartridge.

■ Preparative size

Product number	Particle size (μm)	Inner diameter (mm)	Length (mm)
80214	5	10	150
80215			250
80224		20	150
80225			250
80234		30	150
80235			250

7 References

- High performance liquid chromatography handbook 2nd rev, Maruzen, JSAC Kanto
 Ekikuro-toranomaki, Hiroshi Nakamura et al., Tsukuba press
 Ekikuro-hyounomaki, Hiroshi Nakamura et al., Tsukuba press
 O.Shirota, Y. Ohtsu and O.Nakata, *J.Chromatogr., Sic.*, 28, 553 (1990)

8 Global distributors

South Korea

▶ **Young Jin Biochrom Co., Ltd.**

Kranz Techno Bd. 1301, 5442-1, SangDaeWon-Dong,
JungWon-ku, SungNam-City, 462-120

Tel : 82-2-575-1494

Taiwan

▶ **Analab Corporation**

6F-2, No. 181, Sec. 2, Changan E. Rd., Taipei

Tel : 886-2-27766931

Singapore

▶ **Alpha Analytical (S) PTE LTD**

30 Toh Guan Road #08-02 OSIM Distribution Centre,
Singapore 608840

Tel : (65)6567-8885

Australia

▶ **Phenomenex**

PO Box 4084, Lane Cove, NSW 2066

Tel : 02-9428-6444

USA

▶ **Phenomenex**

411 Madrid Ave. Torrance, CA, 90501

Tel : (310)212-0555

USA

▶ **JM Science**

Grand Island Research Park, PO Box 250/355 Lang
Blvd. Grand Island, NY 14072-0250

Tel : (716)774-8706

UK and Ireland

▶ **Phenomenex**

Melville House, Queens Avenue, Hurdsfield Industrial
Estate, Macclesfield, Cheshire SK10 2BN

Tel : 01625-501367

France

▶ **Interchim**

213 Avenue J.F. Kennedy BP - No 1140

03103 MONTLUCON Cedex

Tel : (33) 470 03 88 55

Belgium, Luxembourg, and Netherlands

▶ **Analis**

Leeuwerikstraat 28, 9000 Gent

Tel : 09/243 77 10

China

▶ **Techmate Ltd.**

C2-F2, 168# FengBao Rd., Feng-Tai DIST, Beijing
100070

Tel : 86-10-8361-3215

India

▶ **Amkette Analytics Ltd.**

B-27, shri Ram Ind. Estate, 13, G.D. Ambekar Road,
Wadala, Mumbai-400 031

Tel : 91-22-2416-1544

Malaysia

▶ **Alpha Analytical (S) PTE LTD**

#47A, Jalan Anggerik Vanilla N 31/N, Kota Kemuning,
40460 Shah Alam, Selangor Darul Ehsan

Tel : (603) 5124 6088

New Zealand

▶ **Phenomenex**

PO Box31-601, Milford 0741, North Shore City

Tel : 09-4780951

USA

▶ **ESA Biosciences**

22 Alpha Road, Chelmsford, MA 01824-4171

Tel : (978)250-7083

USA

▶ **Parma C&S**

25801 Obrero Dr. #1, Mission Viejo CA 92691

Tel : (949)302-8500

Switzerland

▶ **Thermo Fisher Scientific (Schweiz) AG**

Neuhofstrasse 11, CH-4153, Reinach BL

Tel : +41 (61) 716 7700

Germany

▶ **Phenomenex**

Zeppelinstr. 5, 63741 Aschaffenburg

Tel : 06021-58830-0

Italy

▶ **CPS Analytica**

Via Neera, 8/A, 20141 Milano

Tel : (0039) 02 8954-201