


CAPCELL PAK UG (C₁₈, C₈, Ph, CN, NH₂, SCX)

The CAPCELL PAK UG type utilizes a high-purity silica with low metal impurity (<5ppm), that gives a fast separation of basic and polar compounds with sharp symmetrical peaks.

Excellent surface inertness

PYRIDINE/PHENOL TEST


Column : CAPCELL PAK UG120 4.6mm i.d. x 150mm
 Mobile Phase : CH₃CN/H₂O = 30/70
 Flow Rate : 1.0mL/min
 Temperature : 40°C
 Detection : UV 254nm
 Sample : 1. Pyridine
 2. Phenol



C₁₈ UG


Water soluble vitamins

Column : CAPCELL PAK C₁₈ UG120 S5
 4.6mm i.d. x 150mm
 Mobile Phase : (5mmol/L Sodium hexanesulfonate
 + 20mmol/L H₃PO₄, pH2.3) /
 CH₃CN = 91 / 9
 Flow Rate : 1.0mL/min
 Temperature : 40°C
 Detection : UV 210nm
 Sample : 1. L-Ascorbic acid
 2. Nicotinic acid
 3. Nicotinamide
 4. Sodium pantothenate
 5. Pyridoxine hydrochloride
 6. Riboflavin phosphate
 7. Thiamine
 8. Folic acid
 9. Biotin
 10. Riboflavin (VB₂)




Antiseptics

Column : CAPCELL PAK C₁₈ UG120 S5
 4.6mm i.d. x 150mm
 Mobile Phase : 0.05mol/L NaH₂PO₄, pH4.5 / MeOH /
 CH₃CN = 50 / 35 / 15, 4mmol/L
 Cetyltrimethyl ammonium chloride
 Flow Rate : 1.0mL/min
 Temperature : 40°C
 Detection : UV 235nm
 Sample : 1. Methylparaben
 2. p-Hydroxy benzoic acid
 3. Ethylparaben
 4. Dehydroacetic acid
 5. n-Propylparaben
 6. Sorbic acid
 7. Benzoic acid
 8. Iso-Butylparaben
 9. n-Butylparaben
 10. Salicylic acid




Acidic hair dyes

Column	: CAPCELL PAK C ₁₈ UG120 S5 4.6mm i.d. x 150mm
Mobile Phase	: 50mmol/L Sodium 1-octanesulfonate / CH ₃ OH = 60 / 40 (pH2.5, H ₃ PO ₄)
Flow Rate	: 1.0mL/min
Temperature	: 40°C
Detection	: UV 280nm
Sample	: 1. Resorcin 2. p-Aminophenol 3. p-Nitro-o-phenylenediamine 4. m-Aminophenol 5. o-Aminophenol 6. m-Phenylenediamine 7. p-Phenylenediamine 8. p-Toluenediaminosulfonate 9. p-Amino-o-cresol




Ingredients of cold medicines

Column	: CAPCELL PAK C ₁₈ UG120 S5 4.6mm i.d. x 150mm
Mobile Phase	: 0.05mol/L NaH ₂ PO ₄ / CH ₃ CN = 20 / 80 (pH2.5, H ₃ PO ₄)
Flow Rate	: 1.0mL/min
Temperature	: 40°C
Detection	: UV 280nm
Sample	: 1. Potassium guaiacolsulfonate 2. Acetaminophen 3. Caffeine 4. Salicylamide 5. Chlorpheniramine maleate 6. Phenol(I.S.) 7. Aspirin 8. Ethenzamide




Metabolites of testosterone

Column	: CAPCELL PAK C ₁₈ UG120 S3 4.6mm i.d. x 150mm
Mobile Phase	: H ₂ O / CH ₃ OH / THF = 55 / 38 / 7
Flow Rate	: 1.0mL/min
Temperature	: 40°C
Detection	: UV 254nm
Sample	: 1. 6 α -Hydroxytestosterone 2. 7 α -Hydroxytestosterone 3. 6 β -Hydroxytestosterone 4. 16 α -Hydroxytestosterone 5. 16 β -Hydroxytestosterone 6. 2 α -Hydroxytestosterone 7. 2 β -Hydroxytestosterone 8. Androsterone 9. Methyltestosterone




Active ingredients of bearberry leaf

Column	: CAPCELL PAK C ₁₈ UG120 S5 4.6mm i.d. x 150mm
Mobile Phase	: H ₂ O / CH ₃ OH / 0.1mol/L HCl = 94 / 5 / 1
Flow Rate	: 0.6mL/min
Temperature	: 20°C
Detection	: UV 280nm
Sample	: 1. Arbutin 2. Hydroquinone 3. Gallic acid




Tryptic digest of casein (peptide mapping)

Column	: CAPCELL PAK C ₁₈ UG120 S5 1.0 mm i.d. x 250 mm
Mobile Phase	: A : 0.1vol% TFA, H ₂ O B : 0.1vol% TFA, CH ₃ CN B 10% - 50% (45min) Gradient
Flow Rate	: 70 μ L/min
Temperature	: 35°C
Detection	: UV 210nm
Sample	: 1. Tryptic digest of casein



C₈ UG120

Suitable for quick separation of polar compounds, which used to be retained too long in other C₁₈ phases.




Column	: CAPCELL PAK C ₈ UG120 S5 4.6mm i.d. x 150mm
Mobile Phase	: 0.1vol% H ₃ PO ₄ / CH ₃ OH = 70 / 30
Flow Rate	: 1.0mL/min
Temperature	: 40°C
Detection	: UV 254nm
Sample	: 1. Hydroquinone 2. Resorcin 3. Catechol 4. Phenol

Ph UG120

Used for obtaining a different selectivity for analytes possessing an aromatic moiety.

Antiepileptics


Column	: CAPCELL PAK Ph UG120 S5 4.6mm i.d.x 150mm
Mobile Phase	: { 50mmol/L Na ₂ HPO ₄ + 50mmol/L KH ₂ PO ₄ (pH6.8) } / CH ₃ CN = 70 / 30
Flow Rate	: 1.0mL/min
Temperature	: 40°C
Detection	: UV 254nm
Sample	: 1. Phenobarbital 2. Carbamazepine 3. Phenytoin

**CN UG120**

A phase having the least retentive nature of all reversed phases and a different selectivity brought by cyano groups.

Steroids

Column	: CAPCELL PAK CN UG120 S5 4.6mm i.d. x 150mm
Mobile Phase	: CH ₃ CN / H ₂ O = 35 / 65
Flow Rate	: 1.0mL/min
Temperature	: 35°C
Detection	: UV 242nm
Sample	: 1. Cortisol 2. Cortisone 3. Corticosterone 4. Testosterone






NH₂ UG80

To be used as a normal phase under a water/organic mobile phase, or a weak anion exchanger under an acidic buffer.

Nucleotides


Column	: CAPCELL PAK NH ₂ UG80 S5 4.6mm i.d. x 250mm
Mobile Phase	: 0.05mol/L (NH ₄) ₂ HPO ₄ (pH3.0)
Flow Rate	: 1.0mL/min
Temperature	: 40°C
Detection	: UV 254nm
Sample	: 1. 5'-CMP 2. 5'-AMP 3. 5'-UMP 4. 5'-IMP 5. 5'-GMP



SCX UG80

A strong cation exchanger used for basic compounds.


Determination of biogenic amines




Column	: CAPCELL PAK SCX UG80 4.6mm i.d. x 150mm
Mobile Phase	: 0.2mol/L CH ₃ COONa + 0.02mol/L CH ₃ COOH (pH5.6)
Flow Rate	: 1.0mL/min
Temperature	: 40°C
Detection	: UV 254nm
Sample	: 1. Dopa 2. Creatinine 3. Norepinephrine 4. Epinephrine 5. Octopamine 6. Normetanephrine 7. Dopamine 8. Isoproterenol 9. Metanephrine 10. Tyramine 11. Serotonin

Excellent lot-to-lot reproducibility

Column	: CAPCELL PAK SCX UG80 4.6mm i.d. x 150mm
Mobile Phase	: 0.2mol/L NH ₄ H ₂ PO ₄ (pH3.5)
Flow Rate	: 1.0mL/min
Temperature	: 40°C
Detection	: UV 254nm
Sample	: 1. Uracil 2. Thymine 3. Guanine 4. Cytosine 5. Adenine



Highly durable



Column	: CAPCELL PAK SCX UG80 4.6mm i.d. x 150mm
Mobile Phase	: 0.2mol/L NH ₄ H ₂ PO ₄ (pH3.5)
Flow Rate	: 1.0mL/min
Temperature	: 40°C
Detection	: UV 254nm
Sample	: 1. Uracil 2. Thymine 3. Guanine 4. Cytosine 5. Adenine