

CHROMATOGRAPHY IS OUR BUSINESS

# Crimper and Decapper



Easy Grip Manual Crimper & Decapper Electronic Crimper & Decapper High Power Electronic Crimper & Decapper Classic Manual Crimper & Decapper

www.bgb-shop.com

#### BGB Easy Grip Manual Crimper and Decapper

Our Easy Grip Manual Crimper and Decapper for all 8, 11, 13 and 20 mm crimp caps. The curved handle features comfortable use compared to metal grip designs. The bottom-pull handle allows for a steadier hold than top-push handles.

An easily-viewed adjustment knob is located at the head of the crimper and decapper. The clearly shown + and - symbols along with the directional arrows, simplify adjusting the intensity of the crimp or decap needed. The knob is also designed to indicate when the crimp setting has been reached, when it fully contacts the body of the crimper/decapper.



All Easy Grip Manual Crimper and Decapper are equipped with an adjustment knob.

| BGB<br>LAY GRIP<br>charter<br>trans                      |   |   |
|--|---|---|
| Part No.   | Description   | Otv   |
| Part No.   | Description   | Qty.  |
| 295000   | 8 mm Crimper  | ea.   |
|  | •   | -   |
| 295000   | 8 mm Crimper  | ea.   |
| 295000<br>295001   | 8 mm Crimper<br>11 mm Crimper   | ea.<br>ea.                                    |
| 295000<br>295001<br>295002                               | 8mm Crimper<br>11mm Crimper<br>13mm Crimper   | ea.<br>ea.<br>ea.                             |
| 295000<br>295001<br>295002<br>295003                     | 8 mm Crimper<br>11 mm Crimper<br>13 mm Crimper<br>13 mm Crimper for Flip Off Caps   | ea.<br>ea.<br>ea.<br>ea.                      |
| 295000<br>295001<br>295002<br>295003<br>295004           | 8mm Crimper<br>11mm Crimper<br>13mm Crimper<br>13mm Crimper for Flip Off Caps<br>20mm Crimper   | ea.<br>ea.<br>ea.<br>ea.                      |
| 295000<br>295001<br>295002<br>295003<br>295004<br>295005 | 8 mm Crimper<br>11 mm Crimper<br>13 mm Crimper<br>13 mm Crimper for Flip Off Caps<br>20 mm Crimper<br>20 mm Crimper for Flip Off Caps | ea.<br>ea.<br>ea.<br>ea.<br>ea.<br>ea.<br>ea. |



# BGB Classic Manual Crimper and Decapper

- Chemically resistant surface finish especially designed for the application in the labs.
- Hardened crimping jaws made of special alloy that gurantees long life
- Adjustable in crimping pressure limitation by a screw in the handle
- Additionally 11 mm, 13 mm and 20 mm crimpers are adjustable in crimping height by screwing up or down the pressure block in the crimping head with a hexagon key.



| Part No. | Description                                | Qty. |
|----------|--|------|
| 295020   | 8mm Crimper                                | ea.  |
| 295021   | 11 mm Adjustable Crimper                   | ea.  |
| 295022   | 13 mm Adjustable Crimper                   | ea.  |
| 295023   | 13 mm Adjustable Crimper for Flip Off Caps | ea.  |
| 295024   | 20 mm Adjustable Crimper                   | ea.  |
| 295025   | 20 mm Adjustable Crimperf or Flip Off Caps | ea.  |
| 295026   | 8mm Decapper                               | ea.  |
| 295027   | 11 mm Decapper                             | ea.  |
| 295028   | 13 mm Decapper                             | ea.  |
| 295029   | 20 mm Decapper                             | ea.  |

### BGB Electronic Crimper and Decapper Model 6

Our new Electronic Crimpers and Decappers with enhanced OLED screen displays fit on standard laboratory vials and seals with aluminum caps. The electronically adjustable cycle is accurate and reproducible and can be adjusted on the display of the unit. The tools are ergonomically designed to reduce strain and arm injury associated with many manual crimping tools.

- Powerful motor with fast cycle
- Longer lasting battery with Lithium Ion technology
- Advanced charging system including warning light when recharging is necessary
- Can be operated while charging
- Controls on top of unit with enhanced OLED Screen





OLED screen display

| Part No. | Description                             | Qty. |
|----------|---|------|
| 297040   | 8mm Crimper Model 6                     | ea.  |
| 297041   | 11 mm Crimper Model 6                   | ea.  |
| 297042   | 13 mm Crimper Model 6                   | ea.  |
| 297043   | 13 mm Crimper Model 6 for Flip Off Caps | ea.  |
| 297044   | 20 mm Crimper Model 6                   | ea.  |
| 297045   | 20 mm Crimper Model 6 for Flip Off Caps | ea.  |
| 297046   | 11 mm Decapper Model 6                  | ea.  |
| 297047   | 13 mm Decapper Model 6                  | ea.  |
| 297048   | 20 mm Decapper Model 6                  | ea.  |



Note: Not suitable for all-steel or magnetic ring crimp caps, please use High Power Electronic Crimping Tool (see page 4) for this type of caps.

### **BGB** Crimping Tool Stand

Stand for BGB Electronic and Easy Grip Manual Crimpers and Decappers.





Examples of use.

# BGB High Power Electronic Crimping Tool HPS Model 6

The High Power Electronic Crimping Tool with new OLED screen display is faster and more powerful than the electronic crimping tool. The HPS crimping tool is designed to store separate programs for different caps and seals. It is powered by a universal power supply and uses interchangeable jaw sets. Jaw sets are strong enough to crimp all steel and magnetic caps.

- High-Power version of electronic crimping tool with fixed power supply and cord (no battery).
- Uses interchangeable jaw sets.
- Stores separate programs for different caps and seals.
- Strong enough for all steel and magnetic caps.
- Optional base with mounting kit.



#### Jaw Sets for BGB High Power Electronic Crimping Tool HPS



20 mm Crimper for Flip Off Caps 295067 ea.





Jaw Sets are easy to change. Simply insert the jaw set into the bushing at the bottom of the tool. Push up and twist until the set locks into position. To remove, push the outside button and rotate.

DR-BGB-2407

# BGB GC LC MS CE

www.bgb-shop.com

| Switzerland | BGB Analytik AG • Rohrmattstrasse 4 • 4461 Böckten • Phone +41 61 991 00 46 • Fax +41 61 991 00 25 • sales@bgb-analytik.com                         |
|-------------|---|
| Benelux     | BGB Analytik Benelux B.V. • Drielandendreef 42 • 3845 CA Harderwijk • Phone +31 341 700270 • Fax +31 341 700271 • sales.benelux@bgb-analytik.com    |
| France      | BGB Analytik France S.A.S. • 81 Vie de l'Etraz • 01630 St. Jean de Gonville • Phone +33 450 488567 • Fax +33 450 562378 • sales.fr@bgb-analytik.com |
| Germany     | BGB Analytik Vertrieb GmbH • Mühlestraße 1 • 79539 Lörrach • Phone +49 7621 5884270 • Fax +49 7621 5884289 • sales.de@bgb-analytik.com              |
| Turkey      | BGB Analytik Laboratuvar Malzemeleri Ticaret LTD. STİ • Starport Residence • Yenişehir Mah. Osmanlı Bulvarı. • Sümbül Sok. No:10, Daire: 193 •      |
|             | 34912 Kurtköy, Pendik, İstanbul • Phone +90 2169092048 • sales.tr@bgb-analytik.com  |
| USA         | BGB Analytik USA LLC • 8407 F Richmond Hwy • Alexandria, VA 22309 • Phone 1-703-780-1500 • Fax 1-703-991-9131 • sales.usa@bgb-analytik.com          |





11 mm Decapper 295068

ea.

Prices may change without notice and do not include VAT.



295069

ea.

20 mm Decapper 295070 ea.



13 mm

Crimper

295065

ea.

for Flip Off Caps

20 mm

Crimper

295066 ea.